GB Attention Operating Surgeon

device may result in patient injury.

ADVERSE REACTIONS

Complications can potentially occur following any joint replacement surgery. The following complications have occurred in some patients after their hip joint replacement surgery:

- Hip ball (femoral head) and hip socket (acetabulum) may separate (hip dislocation);
- Device loosening from the surrounding bone;
- Allergic reaction to the implant's materials;
- Audible sounds during motion;
- Premature wear or breakage of the implants;
- Bone loss around the implant;
- Periprosthetic fracture;
- Change in the length of the treated leg;
- Hip pain and stiffness;
- Loss of hip flexibility of the hip joint;
- Nerve damage;
- Calcification.

STERILITY

The MiniHip stem implants are supplied sterile. Check the integrity of the packaging carefully. Do not use if the packaging is open or damaged. Do not resterilise. For single use only. The product is not labeled "pyrogen free".

INSTRUMENTS

Instruments used to implant the device are supplied non-sterile and must be sterilised prior to use. Instruments should be cleaned according to the Corin cleaning instructions supplied with the instrument set. After cleaning and prior to sterilisation, the reusable device(s) should be double-wrapped or packaged in FDA cleared (USA only) wraps or pouches. Wrapping should be performed using the appropriate wrapping method (e.g. AAMI CSR wrapping technique). The following sterilisation method has been validated, based on AAMI/ANSI/ISO guidelines and recommendations:

- Method: Moist-Heat Sterilisation
- Cycle: Pre-Vacuum (Pre-Vac)
- Temperature: 270° F (132° C)
- Exposure Time: 4 minutes
- Dry-Time: 30 minutes (minimum, in chamber)
- Cool-Time: 60 minutes (minimum, at room temperature)

In the USA other sterilization methods for instruments should not be used as they have not been validated by Corin.

Failure of the implant may result from: selection of improper components, mal-alignment or mal-positioning of the components, inadequate fixation of components, failure to remove surgical debris prior to closure. The risk of fracture, particularly of smaller sized stems, is higher in patients who are young, physically active, and/

or heavy. Surgeons should give consideration to bone and muscle quality. Clean gloves should be worn when handling implants. Avoid scratching or denting implant surfaces. Do not modify implants. Do not reuse implants, an implant may appear undamaged but previous stress may have created imperfections that would reduce the service life of the implant. Surgeons should be thoroughly familiar with the Corin MiniHip surgical technique, the implants, and the instruments prior to performing the surgery.

Corin designs specialised instruments for their joint replacement systems to aid in the accurate implantation of the prostheses. The use of instruments or implant components from other systems can result in inaccurate fit, sizing, and device failure. Intra-operative fracture or breaking of instruments has been reported. Surgical instruments are subject to wear with normal usage. Instruments which have experienced extensive use or excessive force are susceptible to fracture. Surgical instruments should only be used for their intended purpose.

MR SAFETY INFORMATION

The MiniHip stem has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the MiniHip Stem in the MR environment is unknown. Scanning a patient who has this

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advised to validate whichever method they deem appropriate at their institution. EtO sterilisation and cold sterilisation techniques are not recommended.

The information contained in this package insert was current on the date it was printed but the package insert may have been revised after that date. To obtain a current package insert, please contact Corin.

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