2018-03-15

Medtronic

ENGLISH **IMPORTANT INFORMATION ON THE CLYDESDALE® SPINAL SYSTEM**

PURPOSE

This device is a PEEK (POLYETHERETHERKETONE) interbody fusion device intended for stabilization use and to promote bone fusion during the normal healing process following surgical correction of disorders of the spine. The product should be implanted only by a physician who is thoroughly knowledgeable in the implant's material and surgical aspects and who has been instructed as to its mechanical and material applications and limitations.

DESCRIPTION

The CLYDESDALE® Spinal System consists of PEEK cages of various widths and heights, which include Tantalum markers. These devices can be inserted between two lumbar or lumbosacral vertebral bodies to give support and correction during lumbar interbody fusion surgeries. The hollow geometry of the implants allows them to be packed with autogenous bone graft.

Implied warranties of merchantability and fitness for a particular purpose or use are specifically excluded.

INDICATIONS

The CLYDESDALE® Spinal System is designed to be used with autogenous bone graft to facilitate interbody fusion and is intended for use with supplemental fixation systems cleared for use in the lumbar spine. The CLYDESDALE® Spinal System is used for patients diagnosed with degenerative disc disease (DDD) at one or two contiguous levels from L2 to S1. These DDD patients may also have up to Grade 1 Spondylolisthesis or retrolisthesis at the involved levels. DDD is defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment. **These implants may be implanted via a minimally invasive lateral approach.**

CONTRAINDICATIONS

This device is not intended for cervical spine use.

Contraindications include, but are not limited to:

- Infection local to the operative site
- Signs of local inflammation
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Mental illness
- Any other condition which would preclude the potential benefit of spinal implant surgery, such as the presence of tumors or congenital abnormalities, fracture local to the operating site, elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or a marked left shift in the WBC differential count.
- Suspected or documented allergy or intolerance to composite materials
- Any case not needing a fusion
- Any case not described in the indications
- Any patient unwilling to cooperate with postoperative instructions
- Patients with a known hereditary or acquired bone friability or calcification problem
- Pediatric cases or where the patient still has general skeletal growth
- Spondylolisthesis unable to be reduced to Grade 1
- Any case where the implant components selected for use would be too large or too small to achieve a successful result
- Any case that requires the mixing of metals from two different components or systems
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance
- Prior fusion at the level to be treated

NOTA BENE: Although not absolute contraindications, conditions to be considered as potential factors for not using this device include:

- **Severe bone resorption.**
- **Osteomalacia**
- **Severe osteoporosis.**

POTENTIAL ADVERSE EVENTS

Adverse effects may occur when the device is used either with or without associated instrumentation.

The potential risk of adverse effects as a result of movement and non-stabilization may increase in cases where associated complementary support is not employed. Potential adverse events include but are not limited to:

- Implant migration.
- Breakage of the device(s).
- Foreign body reaction to the implants including possible tumor formation, auto immune disease, and/or scarring.
- Pressure on the surrounding tissues or organs.
- Loss of proper spinal curvature, correction, height, and/or reduction.
- Infection.
- Bone fracture or stress shielding at, above, or below the level of surgery.
- Non-union (or pseudoarthrosis).
- Loss of neurological function, appearance of radiculopathy, dural tears, and/or development of pain.
- Neurovascular compromise including paralysis, temporary or permanent retrograde ejaculation in males, or other types of serious injury.
- Cerebral spinal fluid leakage.
- Hemorrhage of blood vessels and/or hematomas.
- Discitis, arachnoiditis, and/or other types of inflammation.
- Deep venous thrombosis, thrombophlebitis, and/or pulmonary embolus.
- Bone graft donor site complication.
- Inability to resume activities of normal daily living.
- Early or late loosening or movement of the device(s).
- Urinary retention or loss of bladder control or other types of urological system compromise.
- Scar formation possibly causing neurological compromise or compression around nerves and/or pain.
- Fracture, microfracture, resorption, damage, or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above, and/or below the level of surgery.
- Retropulsed graft.
- Herniated nucleus pulposus, disc disruption, or degeneration at, above, or below the level of surgery.
- Loss of or increase in spinal mobility or function.
- Reproductive system compromise, including sterility, loss of consortium, and sexual dysfunction.
- Development of respiratory problems (e.g., pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.)
- Change in mental status.
- Cessation of any potential growth of the operated portion of the spine.
- Death.

WARNINGS AND PRECAUTIONS

A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where other patient conditions may compromise the results. Use of this product without bone graft or in cases that do not develop a union will not be successful.

Preoperative and operating procedures, including knowledge of surgical techniques, good reduction, and correct selection and placement of the implants are important considerations in the successful utilization of the system by the surgeon. Further, the proper selection and the compliance of the patient will greatly affect the results. Patients who smoke have been shown to have a reduced incidence of bone fusion. These patients should be advised of this fact and warned of this consequence. Obese, malnourished, and/or alcohol/drug abuse patients and those with poor muscle and bone quality and/or nerve paralysis are also poor candidates for spinal fusion.

Patients with previous spinal surgery at the levels to be treated may have different clinical outcomes compared to those without a previous surgery.

A device that has been implanted should never be reused, reprocessed, or resterilized under any circumstances. Sterile packaged devices should also never be resterilized. Reuse, reprocessing, or resterilization may compromise the structural integrity of these implants and create a risk of contamination of the implants which could result in patient injury, illness, or death.

PHYSICIAN NOTE: Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

! USA FOR US AUDIENCES ONLY

Caution: Federal law (USA) restricts these devices to sale by or on the order of a physician.

MRI INFORMATION

The CLYDESDALE® Spinal System has not been evaluated for safety, compatibility, heating, or migration in the MR environment.

IMPLANT SELECTION

The selection of the proper size, shape, and design of the implant for each patient is crucial to the success of the procedure. Surgical implants are subject to repeated stresses in use, and their strength is limited by the need to adapt the design to the human anatomy. Unless great care is taken in patient selection, placement of the implant, and postoperative management to minimize stresses on the implant, such stresses may cause material fatigue and consequent breakage or loosening of the device before the fusion process is complete, which may result in further injury or the need to remove the device prematurely.

DEVICE FIXATION

Installation and positional adjustment of implants must only be done with special ancillary instruments and equipment supplied and designated by Medtronic. In the interests of patient safety, it is therefore recommended that Medtronic implants are not used with devices from any other source.

Never, under any circumstances, reuse a CLYDESDALE® Spinal System device. Even when a removed device appears undamaged, it may have small defects or internal stress patterns that may lead to early breakage.

PREOPERATIVE

- Only patients that meet the criteria described in the indications should be selected.
- Patient conditions and/or predispositions such as those addressed in the aforementioned contraindications should be avoided.
- Care should be taken in the handling and storage of the device(s). They should not be scratched or damaged. Devices should be protected during storage especially from corrosive environments.
- Further information about this system will be provided upon request.
- The surgeon should be familiar with the various devices before use and should personally verify all devices are present before the surgery begins.
- The size of device for the case should be determined prior to beginning the surgery. An adequate inventory of implant sizes should be available at the time of surgery, including sizes larger and smaller than those expected to be used.
- Unless supplied sterile, all devices should be cleaned and sterilized before use. Additional sterile components should be available in case of any unexpected need.

INTRAOPERATIVE

- The instructions in any available CLYDESDALE® Spinal System surgical technique manual should be carefully followed.
- At all times, extreme caution should be used around the spinal cord and nerve roots. Damage to the nerves will cause loss of neurological function.
- Breakage, slippage, or misuse of instruments or implants may cause injury to the patient or operative personnel.
- To ensure proper fusion below and around the location of the fusion, autogenous bone graft must be used.
- Bone cement should not be used because this material may make removal of these components difficult or impossible. The heat generated from the curing process may damage or deform the PEEK devices.

POSTOPERATIVE

The physician's postoperative directions and warnings to the patient and the corresponding patient compliance are extremely important.

- Detailed instructions on the use and limitations of the device should be given to the patient. The patient must be warned that loosening, and/or breakage of the device(s) are complications which may occur as result of early or excessive weight-bearing, muscular activity, or sudden jolts or shock to the spine.
- The patient should be advised not to smoke or consume excess alcohol during the period of the bone fusion process.
- The patient should be advised of the inability to bend at the point of spinal fusion and taught to compensate for this permanent physical restriction in body motion.
- It is important that immobilization of union is established and confirmed by roentgenographic examination. If a non-union develops or if the components loosen, migrate, and/or break, the devices should be revised and/or removed immediately before serious injury occurs.
- CLYDESDALE® Spinal System implants are interbody devices and are intended to stabilize the operative area during the fusion process.
- Any retrieved devices should be treated in such a manner that reuse in another surgical procedure is not possible.

PACKAGING

Devices may be supplied in a sterile or non-sterile form. Packages for each of the components should be intact upon receipt. Once the seal on the sterile package has been broken, the product should not be re-sterilized. If a loaner or consignment system is used, all sets should be carefully checked for completeness and all components, including instruments, should be carefully checked to ensure there is no damage prior to use. Damaged packages or products should not be used, and should be returned to Medtronic.

CLEANING

Disassembly/reassembly and cleaning instructions can be found at <http://manuals.medtronic.com/>. Refer to the "Reprocessing Instructions for the Lateral Inserters – M708348B087" for disassembly and cleaning instructions

specific to the DL Inserter instrument (part number 2942001). Refer to the “Important Information for Medtronic Reusable Instruments - 0380035” for cleaning instructions for CLYDESDALE® Spinal System trials.

STERILIZATION

The CLYDESDALE® Spinal System implants are supplied in a sterile form. Never autoclave or in any other way attempt to re-sterilize or re-use the CLYDESDALE® Spinal System implants. Instruments are provided non-sterile and must be sterilized by the hospital prior to surgical use. Remove all packaging materials prior to sterilization. Only sterile products should be placed in the operative field. Unless specified elsewhere, these products are recommended to be steam sterilized by the hospital using one of the following sets of process parameters:

Table 1: Sterilization Cycle Parameters for the United States and Its Territories

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME	MINIMUM DRY TIME ¹
Steam	Gravity Displacement	250°F (121°C)	30 Minutes	30 Minutes
Steam	Gravity Displacement	270°F (132°C)	15 Minutes	30 Minutes
Steam	Gravity Displacement	275°F (135°C)	10 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	270°F (132°C)	4 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	275°F (135°C)	3 Minutes	16 Minutes

For Medical Facilities Located Outside the United States and Its Territories: Some non-U.S. health care authorities recommend sterilization according to these parameters so as to minimize the potential risk of transmission of Creutzfeldt-Jakob disease, especially of surgical instruments that could come into contact with the central nervous system.

Table 2: Sterilization Cycle Parameters for Medical Facilities Outside the United States and Its Territories

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME	MINIMUM DRY TIME ¹
Steam	Gravity Displacement	273°F (134°C)	20 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	273°F (134°C)	4 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	273°F (134°C)	20 Minutes	30 Minutes

¹ The minimum dry times were validated using sterilizers having vacuum drying capabilities. Drying cycles using ambient atmospheric pressure may require longer dry times. Refer to the sterilizer manufacturer’s recommendations.

NOTE: Because of the many variables involved in sterilization, each medical facility should calibrate and verify the sterilization process (e.g., temperatures, exposure times) used for their equipment.

The sterilization cycles listed in Table 2 are not considered by the Food and Drug Administration (FDA) to be standard sterilization cycles. It is the end user’s responsibility to use only sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization cassettes) that have been cleared by the FDA for the selected sterilization cycle specifications (time and temperature).

SERVICING

Inspect all instruments prior to use. Return the instrument to Medtronic if any of the following are observed: corrosion, discoloring, pitting, or any other signs of wear.

Inspect the threaded shaft of the inserter instrument. Return the instrument to Medtronic if threads are damaged or distorted or if the shaft appears bent.

Inspect the silicone handle of the inserter instrument. Return the instrument to Medtronic if the silicone handle is discolored, cut, or damaged in any way.

PRODUCT COMPLAINTS

Any health care professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness, and/or performance, should notify the distributor or Medtronic. Further, if any of the implanted spinal system component(s) ever malfunctions (i.e., does not meet any of its performance specifications or otherwise does not perform as intended), or is suspected of doing so, the distributor should be notified immediately. If any Medtronic product ever malfunctions and may have caused or contributed to the death or serious injury of a patient, the distributor should be notified immediately by telephone, fax, or written correspondence. When filing a complaint, provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint, and notification of whether or not a written report from the distributor is requested.

FURTHER INFORMATION

Recommended directions for use of this system (surgical operative techniques) are available at no charge upon request. If further information is needed or required, contact Medtronic.

Covered by one or more of U.S. Pat. Nos. 8,251,997; 5,860,973; 6,991,654; 7,125,425; and other pending patent applications.

©2018 Medtronic Sofamor Danek USA, Inc. All rights reserved.

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