Dansk	3
Deutsch	7
English	11
Español	15
Français	19
Hrvatski	23
Italiano	27
Nederlands	31
Norsk	35
Polski	39
Suomi	43
Svenska	47

BRUGSANVISNING

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CERAMENT™|G er en injicerbar keramisk knoglesubstitut bestående af calciumsulfat, hydroxyapatit og gentamicinsulfat. CERAMENT™|G leveres som en pasta med 17,5 mg gentamicin/ml.

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INSTRUCTIONS FOR USE

DEVICE DESCRIPTION

CERAMENT™|G is an injectable ceramic bone graft substitute, consisting of Calcium sulfate, Hydroxyapatite and Gentamicin sulfate. CERAMENT™|G delivers 17.5 mg Gentamicin/mL paste.

By combining Calcium sulfate and Hydroxyapatite an optimal balance is achieved between implant resorption rate and bone in growth rate. Calcium sulfate acts as a resorbable carrier for Hydroxyapatite. Hydroxyapatite has a slow resorption rate, high osteoconductivity promoting bone in growth and gives long term structural support to the newly formed bone.

By adding Gentamicin, colonization with Gentamicin sensitive microorganisms can be prevented in order to protect bone healing.

PERFORMANCE

The injecting performance is dependent upon the desired working consistency of the CERAMENT™|G. The optimal injection consistency is obtained at approximately 4 minutes and the injection can be continued until approximately 7 minutes when using a 16 G needle.

Drilling of the material can be performed at 15 minutes.

CERAMENT™|G attains final setting at 20 minutes.

The wound can be closed at approximately 15 minutes.

When using a mold to prepare beads, wait 20 minutes before removing the beads from the mold.

All times are from the start of mixing.

INTENDED USE

CERAMENT™|G is a resorbable ceramic bone graft substitute intended to fill gaps and voids in the skeletal system to promote bone healing.

CERAMENT™|G provides a void/gap filler that during the surgical procedure can augment hardware and bone alignments.

Gentamicin is included in CERAMENT™|G to prevent colonization of Gentamicin sensitive microorganisms in order to protect bone healing.

INDICATIONS

CERAMENT™|G is indicated to be placed into bone voids or gaps in the skeletal system, i.e. extremities, spine, and pelvis, that are not intrinsic to the stability of the bony structure. These osseous defects may be e.g. spontaneously occurring, surgically created, resulting from traumatic injury to the bone, identified during primary surgery and revision surgery, or osseous defects identified around hardware devices.

CONTRAINDICATIONS

- Hypersensitivity to any amino-glycoside antibiotics
- Myasthenia gravis
- Severe renal impairment
- Pre-existing calcium metabolism disorder
- Pregnancy
- Breastfeeding

WARNINGS

- Addition of Gentamicin does not negate the need for systemic antibiotics.

PRECAUTIONS

Practical instruction for sterile / aseptic usage

- Adhere to aseptic technique when preparing Gentamicin solution, please note that the external surface of the vial with mixing liquid is not surface sterilized.
- Adhere to sterile surgical technique when preparing the CERAMENT™|G paste.

Supportive therapy

- Appropriate debridement is needed if pre-existing infection has been identified.
- Control active bleeding and remove blood clots and tissue fragments if open surgery.
- Use of active suction drainage may lead to a decrease in effective dose of Gentamicin.
- Consult and comply with the IFU of any additional utensils used.

Device related

- Contact between CERAMENT™|G and living bone, and provision of normal conditions of fracture healing or bone growth, are prerequisites for good treatment outcome.
- May cause inflammatory reaction if present in soft tissue.
- Not intended for load bearing areas, unless it can be assumed after thorough examination that the cortical bone, surgical fixation, or in situ osteosynthesis is sufficient for load bearing function.
- CERAMENT™|G should not be used in joints since its presence may cause inflammation or mechanical obstruction. When used in the vicinity of joints, care should be taken to ensure that bone voids implanted with CERAMENT™|G are not in continuity with joint spaces, for example by ensuring adequate fracture reduction and/or by using fluoroscopic guidance (as determined by best surgical practice).
- Contact with joint fluid may cause resorption of CERAMENT™|G.
- Overpressurization during injection should be avoided as intra-medullary injection with any bone void filler may lead to fat embolization or embolization of device into the blood stream.
- Do not overfill.
- In Aneurysmal bone cysts (ABCs) and other bone cysts prone to producing large volumes of fluid, there is increased risk of wound drainage, soft-tissue inflammation and wound breakdown if treated by open surgery. Use CERAMENT™|G in bead form rather than complete void filling for these indications.
- No clinical experience with additives in CERAMENT™|G. Using alternative mixing solutions and/or adding other substances to the mixture may affect the product setting in an uncontrolled manner, and may affect the safety and effectiveness of the product.
- When used as an alternative to autograft, please note that like many synthetic bone void fillers CERAMENT™|G is not designed to be osteoinductive or osteogenic.
- If using CERAMENT™|G in conjunction with allograft or autograft, apply each component separately, without intermixing before application, since intermixing may affect the setting time in an uncontrolled manner.
- If using the CERAMENT™|G to augment hardware and bone alignments during surgical procedure, wait until material final setting time for optimal use.

en

Patient related

- Careful examination of patient medical history is recommended.

Related to general use of Gentamicin

- CERAMENT™[G should be used with caution in premature infants because of their renal immaturity, in elderly people and generally in patients with impaired renal function. Diabetes, auditory vestibular dysfunctions, otitis media, a history of otitis media, previous use of ototoxic drugs and a genetically determined high sensitivity to aminoglycoside induced ototoxicity, are other main factors which may pre-dispose the patient to toxicity.
- Concurrent use of other neurotoxic and/or nephrotoxic drugs can increase the possibility of Gentamicin toxicity. Co-administration with the following agents should be avoided:
 - Neuromuscular blocking agents such as succinylcholine and tubocurarine.
 - Other potentially nephrotoxic or ototoxic drugs such as cephalosporins and methicillin.
 - Potent diuretics such as ethacrynic acid and furosemide.
 - Other aminoglycosides.
- To avoid the risk of adverse events, continuous monitoring (before, during and after) of renal function (serum creatinine, creatinine clearance), hepatic and laboratory parameters are recommended.
- Sulfites can cause allergic-type reactions including anaphylactic symptoms and bronchospasm in susceptible people, especially those with a history of asthma or allergy.
- Patient being treated with gentamicin should be under close clinical observation because of its potential toxicity.
- In some patients with impaired renal function, there has been a transient rise in blood-urea-nitrogen, which has usually reverted to normal during or following cessation of therapy.
- Gentamicin should be used with care in conditions characterized by muscular weakness.
- In cases of significant obesity gentamicin serum concentrations should be closely monitored and a reduction in dose should be considered.

SIDE EFFECTS**The following are complications related to the general use of Gentamicin**

- Ototoxicity and nephrotoxicity are the most common side effects associated with Gentamicin therapy. Both effects are related to renal impairment and therefore dosage alteration should be considered in such patients. In addition, there have been rare reports of changes in electrolyte balance including hypocalcaemia and hypokalaemia caused by renal tubular dysfunction. Vestibular damage and ototoxicity is usually reversible.
- As with other aminoglycosides toxicity is related to serum concentration. At serum levels more than 10 µg/mL the vestibular mechanism may be affected and it is advisable to check serum levels to confirm that the initial peak level does not exceed 10 µg/mL and that trough levels do not exceed 2 µg/mL.
- Other adverse reactions associated with Gentamicin therapy include acute renal failure, nausea, vomiting, urticaria, reversible granulocytopenia, hypersensitivity, anaphylactic reactions, anemia, blood dyscrasia, convulsions, central nervous toxicity, abnormal hepatic function, hypomagnesaemia, stomatitis, purpura, allergic contact sensitization and neuromuscular blockade.
- Combinations of antibiotics containing Gentamicin have been associated with rare reports of *Clostridium difficile* diarrhea.

Potential interactions with concomitant use of pharmaceutical agents

- Antibacterials: increased risk of nephrotoxicity with *cephalosporins notably cephalothin*.
- Gentamicin has been known to potentiate anticoagulants such as warfarin and phenindione.
- Antifungals: increased risk of nephrotoxicity with *amphotericin*.
- Cholinergics: antagonism of effect of *neostigmine* and *pyridostigmine*.
- Cyclosporin, cisplatin: increased risk of nephrotoxicity.
- Cytotoxics: increased risk of nephrotoxicity and possible risk of ototoxicity with *cisplatin*.
- Diuretics: increased risk of ototoxicity with *loop diuretics*.
- Muscle relaxants: effect of non-depolarizing muscle relaxants such as *tubocurarine* enhanced. Neuromuscular blockade and respiratory paralysis have been reported from administration of aminoglycosides to patients who have received curare-type muscular relaxants during anaesthesia.
- Indomethacin possibly increases plasma concentrations of gentamicin in neonates.
- Concurrent use of bisphosphonates may increase the risk of hypocalcaemia.
- Concurrent use of Botulinum Toxin and gentamicin may increase the risk of toxicity due to enhanced neuromuscular block.

Local interactions

Consequences of concurrent use of other locally administered drugs are not known.

The following complications have been reported to result from ceramic bone substitutes

- Calcium based bone void fillers may color wound drainage white. This should not be a concern, however be aware of the risk of infection when drainage occurs.
- There have been reports in the literature on idiosyncratic reactions (laryngospasm and tachyarrhythmia) in children up to the age of 15 treated with ceramic bone substitute containing 75-100% Calcium sulfate and 0-25% Calcium phosphate.

DISCLAIMER

- In cases where it is not possible to establish a sufficient wound closure there might be a risk of skin inflammation reaction and/or prolonged wound drainage.
- Bone fracture and wound complications including hematoma, site drainage, infection and other complications are possible side effects of surgery.

COMPONENTS AND COMPOSITIONS

CERAMENT™

CERAMENT™ CMI

Mixing device pre-filled with ceramic bone substitute, a mixture of Calcium sulfate and Hydroxyapatite.

CERAMENT™ JD

Injection device (accuracy of measuring scale \pm 5%).

Valve

To enable connections between CERAMENT™ CMI and syringes.

Tip Extender

Two Tip Extenders in different lengths to be used with CERAMENT™ JD facilitating paste injection.

GENTAMICIN SOLUTION

BONESUPPORT™ SYRINGE

Syringe with rotating adapter to use when preparing the Gentamicin solution. The red mark indicates the correct volume to take.

CERAMENT™ MIXING LIQUID

Sodium Chloride 9 mg/mL liquid, for dissolving the Gentamicin sulfate.

CERAMENT™ GENTAMICIN

Gentamicin; each dose, 1 mL paste, will contain 17.5 mg Gentamicin.

BONESUPPORT™ DP

A ventilated spike to facilitate easy handling when preparing the Gentamicin solution.

ADDITIONAL UTENSILS NEEDED

- Stopwatch

OPTIONAL UTENSILS

- A cannula or needle with a minimum diameter of 16G
- Bead mold

DIRECTIONS FOR USE

- 1 mL paste contains 17.5 mg Gentamicin (provided as Gentamicin sulfate). Adjust volume of paste for sufficient dose of Gentamicin.
- The recommendation for Gentamicin intravenous is 3-6 mg Gentamicin/kg body weight.
- After the paste is in place, avoid any further adjustment.

STEP BY STEP INSTRUCTIONS FOR PREPARATION OF CERAMENT™ G

The box includes two capsules, one with the CERAMENT™ components, and one with the components for Gentamicin solution. Everything inside the pouches and vials is sterile.

Preparation of Gentamicin solution


When handling the components of the gentamicin solution, adhere to aseptic techniques. The content in the vial with mixing liquid is sterile but not the outer surface of the vial.

1. Open the vial containing CERAMENT™ MIXING LIQUID. Attach the BONESUPPORT™ SYRINGE and withdraw liquid up to the red mark on the syringe.
2. Disconnect the BONESUPPORT™ SYRINGE from the CERAMENT™ MIXING LIQUID vial by separating the rotating adapter and the CERAMENT™ MIXING LIQUID vial from the BONESUPPORT™ SYRINGE.
3. Attach the BONESUPPORT™ DP to the BONESUPPORT™ SYRINGE. Push the syringe with the spike through the membrane of the vial and inject the mixing liquid into the vial and dissolve the Gentamicin powder.
4. Withdraw the Gentamicin solution back into the BONESUPPORT™ SYRINGE. The Gentamicin solution is now ready to be injected into CERAMENT™ CMI.

Preparation of CERAMENT™ G paste

When handling the CERAMENT™ which is supplied in a blister in a pouch, adhere to sterile surgical techniques.

1. Remove the plug on CERAMENT™ CMI.
2. Attach the Valve with the clear end to CERAMENT™ CMI by turning it clockwise and attach the BONESUPPORT™ SYRINGE with Gentamicin solution by pressing the syringe through the blue membrane and turning it clockwise.
3. Inject the Gentamicin solution into the CERAMENT™ CMI. Empty the syringe completely and avoid back flush by keeping the plunger pushed to the bottom before detaching the syringe.
4. Detach the syringe. Retract the blue handle on the CERAMENT™ CMI and remove the red plunger stop.
5. Start the stop watch ($t = 0$ seconds). Mix in a tipped down position for 30 seconds with a frequency of approximately one complete stroke per second. Rotate the blue handle at the end positions.
6. Finish the mixing by:
 - Fully retracting the blue handle into its back position.
 - Lock the plunger by turning the blue collar clockwise 180° until a "click" is heard.
7. Attach the CERAMENT™ JD to the Valve by pressing it through the blue membrane and turning clockwise. Transfer the paste from the CERAMENT™ CMI immediately, with the grading on the CERAMENT™ JD facing the user. When the CERAMENT™ JD is completely filled, excess paste will begin to ooze from under the sleeve. Stop filling when this occurs.
8. Detach the filled CERAMENT™ JD from the Valve and remove the red plunger stop rod. The paste is now ready to use. If applicable, attach Tip Extender or an optional sterile needle (minimum 16G) to the CERAMENT™ JD.

RESTRICTIONS

CERAMENT™|G may only be sold, distributed, and used in accordance with the intended use.

Only for professional use.

The contents of this document may not be duplicated without written permission from BONESUPPORT AB.

CERAMENT™ products are protected by patents:
<http://www.bonesupportpatents.com/>

All trademarks herein are the property of BONESUPPORT AB, or its subsidiaries unless otherwise indicated. Registered trademarks of BONESUPPORT AB:
<https://www.bonesupport.com/trademarks.html>

Store CERAMENT™|G unopened in a clean and dry environment in room temperature (15–25°C / 59–77°F).

Do not use if any of the packages are open or damaged or if the expiration date has been exceeded.

Excess material and opened but unused items must be discarded. Used material should be discarded in accordance with hospital procedures.

STERILITY

CERAMENT™|G is supplied sterile, except for CERAMENT™|MIXING LIQUID that has a non-sterile outer surface. Method for sterilization is by gamma irradiation, steam and ethylene oxide.


CERAMENT™|G is disposable and intended only for single use; the product should not be re-sterilized by any method and should not be re-used due to contamination risks.

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Dokumentets innehåll får inte kopieras utan skriftligt medgivande från BONESUPPORT AB.