BIOMET

Biomet® Modular Femoral Heads
ATTENTION OPERATING SURGEON

5401000257

F

BIOMET UK LIMITED
WATERTON INDUSTRIAL ESTATE
BRIDGEND, CF31 3XA
UNITED KINGDOM
TEL: +44 (0) 1656 655221
FAX: +44 (0) 1656 645454

REF: 5401000257
REVISION: F
DATE: November 2020
CHANGE CONTROL NO: 231223

For additional information please see www.biomet.com

EN

BIOMET® MODULAR FEMORAL HEADS ATTENTION OPERATING SURGEON

IMPORTANT INFORMATION FOR THE OPERATING SURGEON

Before using a product placed on the market by Zimmer or Biomet, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g., product literature, surgical technique). Zimmer or Biomet is not liable for complications that may arise from the use of the device in circumstances outside of Zimmer or Biomet's control including, but not limited to, product selection and deviations from the device's indicated uses or surgical technique.

DESCRIPTION

Biomet modular femoral heads are available with type 1 and type 12/14 tapers with various neck lengths and are compatible with any of Biomet's modular femoral stems of the correct taper type. The modular femoral heads are utilised in conjunction with other hip prosthesis components to form part of a total joint replacement system. Total hip joint prostheses may include: acetabular shells, liners, inserts and cups, femoral stems and femoral heads.

The Exceed ABT M2a metal modular heads are indicated to be used with Biomet's range of ultra high molecular weight acetabular liners and cups.

Metal (CoCrMo) modular femoral heads are indicated for use with Biomet's range of ultra high molecular weight polyethylene acetabular liners and cups.

Biomet has evaluated the compatibility of Biomet devices with implants and components from Zimmer orthopaedic companies. Only authorized combinations should be used. To determine whether these devices have been authorized for use in a proposed combination with Zimmer products, please contact your sales representative and reference the following website: www.productcompatibility.zimmer.com.

MATERIALS

Metal Modular Heads	Cobalt-chromium-molybdenum alloy
Exceed™ ABT™ M2a™ Modular Metal Heads	High carbon cobalt-chromium-molybdenum alloy

INDICATIONS

The indications for the use of hip prosthesis include:

1. Non-inflammatory degenerative joint diseases including osteoarthritis, avascular necrosis and post-traumatic arthritis.
2. Rheumatoid arthritis.
3. Correction of functional deformity.
4. Revision of failed joint reconstruction or treatment.
5. Treatment of femoral neck and trochanteric fractures of the proximal femur with femoral head involvement and which is unmanageable using other techniques.

CONTRAINDICATIONS

Absolute contraindications: infection, sepsis and osteomyelitis.

Relative contraindications:

1. Uncooperative patient or patient with neurologic disorders who are incapable of following instructions.
2. Osteoporosis.
3. Metabolic disorder which may impair bone formation.
4. Osteomalacia.
5. Local and distant foci of infection.
6. Rapid joint destruction, marked bone loss or bone resorption apparent on roentgenogram.
7. Vascular insufficiency, muscular atrophy, or neuromuscular disease.
8. Skeletal immaturity.
9. Morbid obesity.
10. Foreign body sensitivity. Where suspected or unknown, material sensitivity tests are to be made prior to implantation.
11. Any condition that may interfere with the survival of the implants such as Paget's disease, Charcot's disease, sickle cell anaemia or traits, lower extremity muscular atrophy or neuromuscular disease.

PATIENT SELECTION FACTORS

Patient selection factors which should be considered include:

1. Primary need to obtain pain relief.
2. Primary need to improve function.
3. Ability of the patient to follow instructions, including control of weight and activity level, and a willingness to undergo rehabilitation.
4. A good nutritional state of the patient.
5. Whether the patient has reached full skeletal maturity.

WARNINGS AND PRECAUTIONS

1. Clinical outcome may be affected by component positioning. Proper placement of the implant should take into consideration individual patient anatomy as well as surgeon preference. The surgical technique sets forth suggested guidelines for placement of the implant including inclination and anteversion.
2. Improper selection, placement, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may lead to subsequent reduction in the service life of the prosthetic components.
3. Malalignment of components or inaccurate implantation can lead to excessive wear and/or failure of the implant or procedure.

4. Complete preclosure cleaning and removal of bone cement debris, metallic debris and other surgical debris at the implant site is critical to minimise wear of the implant articular surfaces.
5. Improper preoperative or intra-operative implant handling or damage (scratches, dents, etc.) can lead to crevice corrosion, fretting, fatigue fracture, and/or excessive wear.
6. Do not modify implants.
7. The surgeon is to be thoroughly familiar with the implants, instruments and surgical procedure prior to performing surgery.
8. Biomet joint replacement prostheses provide the surgeon with a means of reducing pain and restoring function for many patients. While these devices are generally successful in attaining these goals, they cannot be expected to withstand the same activity levels and loads of normal healthy bone and joint tissue.
9. Specialised instruments are designed for specific Biomet hip systems to assist in the accurate implantation of the implants. Use only with compatible instrumentation to minimize risk of inaccurate fit, sizing, and/or excessive wear and device failure.
10. Whilst rare, intraoperative fracture or breakage of instruments has been reported. Instruments that have experienced excessive use or excessive force are susceptible to fracture. Biomet recommends that instruments are examined for wear or damage prior to surgery.
11. All trial, packaging, and instrument components must be removed prior to closing the surgical site. Do not implant.
12. Do not reuse implantable devices. While an implant/device may appear undamaged, previous stress may have created imperfections that would reduce the service life of the implant/device. Do not treat patients with implants / devices that have been even momentarily placed in a different patient. Furthermore, re-using an implant could cause patient contamination.
13. Use clean surgical gloves when handling implants.
14. Do not use any component that is nicked, scratched, chipped or otherwise altered.
15. Do not use Biomet prosthetic implants with femoral stems or acetabular components offered by other manufacturers. Mismatching of components or taper sizes can be expected to cause intraoperative or postoperative fracture of femoral heads or component dissociation. All modular components must be accurately seated to reduce the risk of dissociation. For specific Zimmer or Biomet product compatibility, reference the following website: www.productcompatibility.zimmer.com
16. Use Biomet femoral and modular head component with appropriate matching "Type I Taper" or "12/14 Taper".
17. Firmly seat modular head components to reduce the risk of dissociation and/or micro motion. Micro motion may contribute to difficult head removal.
18. Thoroughly clean and dry taper prior to attachment of the modular head component to minimise the risk of crevice corrosion and improper seating.
19. Ceramic inserts are not suitable for use with metal femoral heads.
20. Patients should be warned of the impact of excessive loading that can result if the patient is involved in an occupation that includes substantial walking, running, lifting, or excessive muscle loading due to weight that place extreme demands on the hip and can result in device failure or dislocation.
21. Patient smoking may result in delayed healing, non-healing and/or compromised stability in or around the placement site.
22. The appropriate surgical technique will be determined by compatible components used to complete the system. Deviations from the recommended technique could result in early loosening/failure of the device, or other adverse events as outlined in the following section- 'POSSIBLE ADVERSE EFFECTS'. Corrosion is rare but has been reported with the use of dissimilar metals in total hip replacement.

Please refer to the suggested operative technique guide available from your sales representative for suggested operative technique.

MRI COMPATIBILITY

Biomet Hip System in the Magnetic Resonance (MR) Environment

The risks associated with a passive implant in an MR environment have been evaluated and are known to include heating, migration, and image artifacts at or near the implant site.

Non-clinical testing has demonstrated that Zimmer Biomet Total Hip Implant Systems are MR Conditional as indicated by the MR symbol. A patient with this device can be safely scanned in a MR system meeting the following conditions:

MR Information

Safety information for the use of MRI procedures (i.e. imaging, angiography, functional imaging, spectroscopy, etc.) pertains to shielded MRI systems under the following specifications:

- Static magnetic field of 1.5-Tesla (1.5 T) and 3.0-Tesla (3.0 T)
- Maximum spatial gradient field of 1300 Gauss/cm or less when used with a stainless steel hip component and 2500 Gauss/cm or less when used with a cobalt-chromium alloy or titanium alloy hip component.
- Maximum MR System reported, whole-body-averaged specific absorption rate (SAR) of:
 - 2 W/kg for 15 minutes of scanning for patient landmarks above the umbilicus and
 - 1 W/kg for 15 minutes of scanning for patient landmarks below the umbilicus.
- Padding for protection against Radio Frequency (RF) burns should be placed between the wall of the bore and extremities.
- Insulating padding between the knees to prevent legs from touching.
- Arms and hands of the patient should not touch each other or other bare skin.
- Quadrature Transmit Mode Only

The effects of MRI procedures using MR systems and conditions above these levels have not been determined. The health state of the patient or the presence of other implants may require a lowering of MR limits.

MR Heating

Under the scan conditions defined above, the Biomet hip system implants are expected to produce a maximum temperature rise of less than 3°C after 15 minutes of continuous scanning.

Image Artifacts

In non-clinical testing, the image artifact caused by the device extends up to 100 mm from stainless steel components and 80 mm from cobalt-chromium and titanium implants when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

Other

In non-clinical 3.0 T testing, the materials used in Biomet implant systems did not produce any magnetically induced displacement force or torque that would result in migration of the devices in the spatial gradient and static fields identified above.

PATIENT WARNINGS

Postoperative care is important. The patient must be warned of the limitations of the reconstruction and the need for protection of the implants from full weight bearing or load bearing until adequate fixation and healing has occurred. The patient is to be cautioned to govern activities, protecting the joint replacement from unreasonable stress conditions. Excessive activity, failure to control body weight and trauma affecting the joint replacement have been associated with premature failure of the reconstruction by loosening, fracture and/or wear of the implants. Loosening of the components can result in increased production of wear particles, as well as accelerating damage to bone, making successful revision surgery more difficult. The patient is to be advised of the importance of postoperative follow-up examination and that excessive activity or trauma can lead to failure. The patient is to be warned in advance of the surgical risks and made aware of possible adverse effects. The patient must be warned that the device does not replace normal healthy bone, and that the implant can break or be damaged as a result of excessive load bearing or trauma. The patient must be warned to inform any other medical practitioner who may treat him in the future of the presence of the implant.

POSSIBLE ADVERSE EFFECTS

1. Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions involving various sizes of macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to or during the healing process. Particulate wear debris and discoloration from metallic and polyethylene components of joint implants may be present in adjacent tissue or fluid. It has been reported that wear debris may initiate a cellular response resulting in osteolysis, which may result in loosening of the implant.
2. Early or late postoperative infection and allergic reaction.
3. Intraoperative bone perforation or fracture may occur, particularly in the presence of poor bone stock caused by osteoporosis, bone defects from previous surgery, bone resorption, or while inserting the device.
4. Loosening or migration of the implants can occur due to loss of fixation, trauma, malalignment, bone resorption, and/or excessive activity.
5. Periarticular calcification or ossification with or without impediment of joint mobility.
6. Inadequate range of motion due to improper selection or positioning of components.
7. Leg length discrepancy.
8. Dislocation and subluxation due to inadequate fixation and improper positioning. Muscle and fibrous tissue laxity can also contribute to these conditions.
9. Fatigue fracture of component can occur as a result of loss of fixation, strenuous activity, malalignment, trauma, non-union, and/or excessive weight.
10. Fretting and crevice corrosion can occur at interfaces between components.
11. Wear and/or deformation of articulating surfaces.
12. Trochanteric avulsion or non-union as a result of excess muscular tension, early weight bearing, or inadequate reattachment.
13. Problems of the knee or ankle of the affected limb or contralateral limb aggravated by leg length discrepancy, too much femoral medialisation or muscle deficiencies.
14. Intraoperative or postoperative bone fracture and/or postoperative pain.
15. Peripheral neuropathies have been reported following total joint replacement surgery. Subclinical nerve damage occurs more frequently, possibly as a result of surgical trauma.
16. Cardiovascular reactions: venous thrombosis, pulmonary embolism or myocardial infarction.
17. Haematoma.
18. Delayed wound healing.
19. Altered gait/limp.
20. Bursitis.
21. Non-malignant, non-infective soft tissue masses, sometimes referred to as pseudotumours.
22. Aseptic lymphocyte-dominated vasculitis associated lesion (ALVAL).
23. Gastrointestinal or genitourinary complications.
24. Pulmonary embolism.
25. Death.
26. Effusion.
27. Wear and corrosion of the metal components can cause an "adverse local tissue reaction (ALTR)" or an "adverse reaction to metal debris (ARMD)", which can damage the surrounding bone and soft tissues. The debris can cause soft tissue necrosis and adversely affect the results of revision surgery.
28. Pain, swelling, or the onset of a limp may occur, requiring further evaluation from an orthopaedic surgeon.

STERILITY

Unless otherwise specified Biomet joint replacement components are sterilised by gamma irradiation (25 – 40 kGy). This is denoted on the label by the symbol:

STERILE **R**

The packages should be inspected for damage or flaws in the sterile barrier before opening. Product should be assumed to be non-sterile if the packaging is damaged. Do not use if packaging is damaged.

All trademarks herein are the property of Biomet, Inc. or one of its affiliates unless otherwise indicated.

CE Mark on the package insert (IFU) is not valid unless there is a CE Mark on the product (description) label.

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