Bar code
87-6203-792-99

Rx only

Rev. C January 2016
© 2016, 2012, Zimmer, Inc.

**FEMORAL HEADS FOR TAPERED-NECK
FEMORAL HIP STEMS
UNIPOLAR ENDOPROSTHESIS FEMORAL HEADS**

zimmer

Zimmer, Inc.
1800 West Center Street
Warsaw, Indiana 46580
USA
1-800-348-2759 (US only)
+1-574-372-4999

EC REP

Zimmer U.K. Ltd.
9 Lancaster Place
South Marston Park
Swindon, SN3 4FP, UK

ENGLISH

The CE mark is valid only if it is also printed on the product label.

**FEMORAL HEADS FOR TAPERED-NECK FEMORAL HIP STEMS
UNIPOLAR ENDOPROSTHESIS FEMORAL HEADS**

Before using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g., product literature, written surgical technique). Zimmer is not liable for complications that may arise

from the use of the device outside of Zimmer's control including, but not limited to, product selection and deviations from the device's indicated uses or surgical technique.

Due to the acquisition of pre-existing product lines, Zimmer has initiated a testing program to evaluate the compatibility of these devices with implants and components made or distributed by all Zimmer Companies, which includes Zimmer GmbH (previously Centerpulse Orthopedics, Ltd.), Zimmer Trabecular Metal Technology, Inc. (previously Implex Corp.), Zimmer U.K. Ltd., and Zimmer Austin, Inc. (previously Centerpulse Orthopedics, Inc.). Only authorized combinations should be used. To determine whether these devices have been authorized for use in a proposed combination, please contact your Zimmer sales representative or visit the Zimmer website: www.productcompatibility.zimmer.com. A printout of the website information can also be obtained by calling Zimmer, Inc. Customer Service, 1-800-348-2759 (U.S.) or the local international access code +1-574-372-4999 (outside the U.S.).

DESCRIPTION

Metal femoral heads for tapered-neck femoral hip stems are offered in two material types, where available: *Zimaloy*[®] Cobalt-Chromium-Molybdenum Alloy and *Zimtron*[™] High Strength Stainless Steel. All femoral heads are supplied in a variety of diameters and neck lengths to use in total hip arthroplasty.

- Femoral heads manufactured from *Zimaloy* Alloy are intended for mating with *Tivanium*[®] Ti-6Al-4V Alloy and *Zimaloy* Alloy femoral stems equipped with tapered necks. These femoral heads are offered in two Morse-type taper options, 6-degree and 12/14, to mate with the corresponding 6-degree or 12/14 taper of a femoral stem. Additional femoral heads that have been treated with a surface hardening process are also available for specified sizes.
- Femoral heads manufactured from *Zimtron* Stainless Steel are intended for mating with stainless steel femoral stems equipped with tapered necks, where available. These femoral heads are available in two Morse-type taper options, 6-degree and 12/14, to mate with the corresponding 6-degree or 12/14 taper of a femoral stem.

The Unipolar Endoprosthesis Femoral Head is made from *Zimaloy* Alloy and is intended for mating with *Tivanium* or *Zimaloy* Alloy femoral stems equipped with tapered necks. Designs are available for mating with both 6-degree and 12/14 proximal neck tapers.

- The Unipolar Endoprosthesis Femoral Head is available in diameters from 38 to 63mm.
- If a longer neck configuration is needed, the Unipolar Endoprosthesis Head Adapter may be used. The Head Adapter is available in *Tivanium* Alloy for the 6-degree stem taper and *Zimaloy* Alloy for the 12/14 stem taper.

INDICATIONS

Femoral heads are intended for use in total hip arthroplasty; the Unipolar Endoprosthesis Femoral Heads are intended for use in hemi-hip arthroplasty. Refer to the femoral stem package insert for the full description of indications for the construct.

INFORMATION FOR USE

Special instruments are available for securing the head on the femoral stem and for removing the head if its diameter or neck length is not appropriate.

WARNINGS

- Do not use another manufacturer's femoral component with a Zimmer Femoral Head or Unipolar Endoprosthesis Femoral Head/Adapter. These devices are designed and intended to be used only with Zimmer femoral stems.
- The 6-degree and 12/14 neck taper components are not interchangeable with each other.
- Do not couple the Head Adapter with femoral heads designed for total hips.
- Do not attempt removal of the Head Adapter from the Unipolar Endoprosthesis Femoral Head once it has been seated onto the head. Damage to the head, adapter, or both may occur.
- Do not couple a stainless steel component with either cobalt-chrome or titanium alloy components. Galvanic corrosion may result.
- Do not use Femoral or Unipolar Heads for trial reductions. Provisionals are available for this purpose.

- Do not attempt removal of the head from the tapered neck of the femoral stem with any instrument other than the specifically-designed distraction instruments.
- If it becomes necessary to remove the Unipolar Endoprosthesis Femoral Head from the tapered neck of the femoral stem, exercise great care to avoid scratching, notching, or otherwise damaging the Unipolar Head or femoral stem.
- This device is for single use only. Do not reuse.

PRECAUTIONS

Prior to the impaction of the femoral head or unipolar head onto a stem, remove any blood or debris from the Morse-type taper neck and examine the inside of the head to ensure that no particulate matter is present. Failure to ensure that mating surfaces are clean could result in inadequate seating of one component upon the other and subsequent disassembly of the femoral head from the femoral stem.

MAGNETIC RESONANCE (MR) SAFETY AND COMPATIBILITY

In the event of exposure to foreseeable environmental conditions, such as magnetic fields, the user and/or patient should be informed of the following precautions:

- This device has not been evaluated for safety and compatibility in the MR environment.
- This device has not been tested for heating or migration in the MR environment.
- There is a potential for heating and migration in the MR environment.
- There is the potential for metal implants to create MR imaging artifacts in the vicinity of the implant.

ADVERSE EFFECTS

The following adverse effects have been reported:

- Peripheral neuropathies
- Deep wound infections
- Perforation of the acetabulum or femur
- Wear
- Heterotopic bone formation
- Metal sensitivity
- Inflammatory reactions and osteolysis
- Dislocation and subluxation
- Vascular complications
- Trochanteric problems
- Subclinical nerve damage
- Corrosion of metal implants
- Early or late loosening of components
- Pelvic, femoral, or acetabular fractures
- Disassembly of modular components
- Fatigue fracture

STERILITY

Gamma irradiation is indicated by the **STERILE | R** symbol on the labeling. These devices remain sterile as long as the package integrity has not been violated. Inspect each package prior to use and do not use the component if any seal or cavity is damaged or breached or if the expiration date has been exceeded. Once opened, the component must be used, discarded, or resterilized.

STERILIZATION INSTRUCTIONS

- These sterilization instructions are consistent with ANSI/AAMI/ISO standards and guidelines. They should be used for sterile items that were opened but unused.
- Do not reuse instruments or devices labeled for single use only.
- In the event of inadvertent loss of sterility while preparing for surgery, sterile packaged, all metal implants may be resterilized only once for immediate use.. This is subject to the exceptions listed below.
- **DO NOT RESTERILIZE:**
 - Single use only components that have been contaminated with body fluids or debris or previously implanted.
 - Components with a packaging expiration date that has been exceeded.

- Do not use the original plastic cavities or lids for resterilization. Single devices should be packaged in a medical grade sterilization pouch or wrap which conforms to the recommended specifications for steam sterilization provided in the table below. Ensure that the pouch or wrap is large enough to contain the devices without stressing the seals or tearing.
- Aggressive cleaning with detergents and brushes may damage special features of the implant, such as fiber metal pads or bead coatings. Also, certain detergents may be difficult to rinse off polymer items, especially those made of silicone rubber.
- Leave the knitted femoral head cover on during resterilization. Remove it just prior to implantation.
- Modular implant components must be sterilized separately to minimize potential bioburden buildup in the dead space and expansion/contraction stresses.

Recommended Sterilization/Resterilization Specifications

Follow the sterilizer manufacturer's instructions for loading patterns and selection of sterilization parameters. Drying times vary according to load size and should be increased for larger loads.

Single, All Metal Implants

Steam Sterilization

Type	Temperature	Exposure Time	Minimum Dry Time
Gravity Displacement	121°C (250°F)	30 minutes	15 minutes
Gravity Displacement	132°C (270°F)	15 minutes	
UK Pre-vacuum/ Pulsating Vacuum ¹	134°C (273°F)	3 minutes	
Pre-vacuum/Pulsating Vacuum	132°C (270°F)	4 minutes	

¹This cycle is not for use in the United States.

Please contact Zimmer at the following number if you have additional questions. In the USA, call 1-800-348-2759. For calls outside the USA, call the local international access code +1-574-267-6131.

PATIENT COUNSELING INFORMATION

Complications and/or failure of prosthetic implants are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Physical activity or trauma can result in loosening, wear, and/or fracture of the implant. The patient must be counseled about the capabilities of the implant and the impact it will have on his or her lifestyle. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant may not last the rest of the patient's life, or any particular length of time. Because prosthetic implants are not as strong, reliable, or durable as natural, healthy tissues/bones, all such devices may need to be replaced at some point.

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