Barcode 87-6203-750-99	<b>Rx</b> only		Rev. B January 2016 © 2016, 2012, Zimmer, Inc.
<b>TRABECULAR METAL™ ACETABULAR REVISION SYSTEM CAGE</b>			
 <b>zimmer</b>	 Zimmer, Inc. 1800 West Center Street Warsaw, Indiana 46580 USA 1-800-348-2759 (US only) +1-574-372-4999	 Zimmer, UK Ltd. 9 Lancaster Place South Marston Park Swindon, SN3 4FP, UK	 

## ENGLISH

The CE mark is valid only if it is also printed on the product label.

## TRABECULAR METAL™ ACETABULAR REVISION SYSTEM CAGE

Before using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g., product literature, written surgical technique). Zimmer is not liable for complications arising from the use of the device outside of its indicated uses, surgical technique or judgment, product selection, and similar matters outside the control of Zimmer.

- **Zimmer has not tested the safety or effectiveness of this device for use in combination with non-Zimmer products or components. If surgeons elect to assemble and implant a construct that includes components not manufactured or distributed by Zimmer, they do so in reliance on their own clinical judgment and should so inform their patients.**
- **Due to the acquisition of pre-existing product lines, Zimmer has initiated a testing program to evaluate the compatibility of this device with implants and components made or distributed by all Zimmer orthopaedic companies, which include Zimmer GmbH (previously Centerpulse Orthopedics Ltd.), Zimmer, Inc., Zimmer Trabecular Metal Technology, Inc. (previously Implex Corp.), Zimmer U.K. Ltd., and Zimmer Austin, Inc. (previously Centerpulse Orthopedics, Inc.). Only authorized combinations should be used. To determine whether this device has been authorized for use in a proposed combination, please contact your Zimmer sales representative or visit the Zimmer website: [www.productcompatibility.zimmer.com](http://www.productcompatibility.zimmer.com). A printout of the website information can also be obtained by calling Zimmer, Inc. Customer Service, 1-800-348-2759 (U.S.) or the local international access code +1-574-372-4999 (outside the U.S.). Former Centerpulse and Implex products that are now packaged in Zimmer boxes, and for which compatibility could be an issue, have been labeled “former Centerpulse” and “former Implex” to provide clarification for the user.**

## **DESCRIPTION**

The *Trabecular Metal* Acetabular Revision System Cage is a metallic dome-shaped, flanged acetabular component with multiple screw holes for adjunct intraoperative peripheral stabilization of revision shell/augment constructs. The cage is designed to be layered between a *Trabecular Metal* acetabular revision shell and a polyethylene acetabular liner using cement fixation. Additional optional screw fixation is available for the cage and shell. This device is intended to aid in initial stabilization of the *Trabecular Metal* shell until biological fixation can occur. The cage is available in right and left versions and with short and long flanges to meet various anatomical needs. The cage material is Commercially Pure (CP) Titanium.

## INDICATIONS

This device is indicated for patients with conditions of, but not limited to:

- Acetabular dysplasia, osteoporosis, protrusio acetabuli, cystic acetabular roof, reconstruction in cases of defects after fracture, acetabular loosening, tumors or revision surgery.
- Advanced joint destruction resulting from degenerative, post-traumatic, or rheumatoid arthritis.
- Failed previous surgery, e.g., osteosynthesis, joint reconstruction, arthrodesis, hemi-arthroplasty, or total hip replacement.

## INFORMATION FOR USE

- Use only instruments specifically designed for use with this device to help ensure accurate surgical implantation, soft tissue balancing, and evaluation of hip function.

## CONTRAINDICATIONS

- Contraindications include:
  - Osteoradionecrosis
  - Skeletal immaturity
  - Any neuromuscular disease in the affected limb
  - Systemic or local infection

## WARNINGS

- **DO NOT REVERSE BEND FLANGES OF THE *TRABECULAR METAL ACETABULAR REVISION SYSTEM CAGE*.**
- For safe and effective use of this implant, the Surgical Technique for the device should be consulted and carefully followed.
- This device is for single patient use only. Never reuse an implant. Although it may appear undamaged, previous stresses may have created imperfections that could reduce service life of the implant.
- Do not use any component if damage is found or caused during setup or insertion.

- Where there is loss of, or insufficient, acetabular bone stock, bone grafting or other adjunctive reinforcement procedures to provide socket support and cup containment are recommended.
- Stems with protruding features (i.e., collared stems, stems with trochanteric buildup) can increase the potential of metal-to-metal impingement and should be avoided if other alternatives are available.
- Patient physical activity can result in post-operative loosening, wear, and/or fracture of the hip implant.

## **PRECAUTIONS**

- Continued surveillance for new or recurrent sources of infection should be continued as long as the device is in place.
- Do not use implant components to perform a trial reduction of the hip. Do not mate provisionals with implants.
- Ensure that the screw heads are properly seated. Use a 3.2mm diameter drill prior to insertion of the compatible 6.5mm diameter screws (6624, 6250, 03147, 42.19, and 00118 series). Avoid penetration beyond the inner cortex of the pelvis when drilling screw holes and inserting screws.
- Placing screws in the posterior inferior and posterior superior quadrants of the acetabulum may reduce the risk of injury to vascular or neurological structures. A wall chart illustrating these zones is available through your Zimmer representative.
- Appropriate placement of the component helps to avoid unintended loading of the component and may improve service life of the construct.

## **MAGNETIC RESONANCE (MR) SAFETY AND COMPATIBILITY**

In the event of exposure to foreseeable environmental conditions, such as magnetic fields, the user and/or patient should be informed of the following precautions:

- This device has not been evaluated for safety and compatibility in the MR environment.
- This device has not been tested for heating or migration in the MR environment.
- There is a potential for heating and migration in the MR environment.
- There is the potential for metal implants to create MR imaging artifacts in the vicinity of the implant.

## **ADVERSE EFFECTS**

- The following adverse effects have been reported:
  - Peripheral neuropathies
  - Deep wound infections
  - Perforation of the acetabulum or femur
  - Heterotopic bone formation
  - Metal sensitivity
  - Wear
  - Inflammatory reactions and osteolysis
  - Fatigue fracture
  - Vascular complications
  - Trochanteric problems
  - Subclinical nerve damage
  - Corrosion of metal implants
  - Early or late loosening of components of total hip arthroplasty
  - Pelvic, femoral, or acetabular fractures
  - Disassembly of modular components

## **STERILITY**

This device is provided sterile by prior exposure to gamma irradiation (indicated by the

**STERILE R** symbol on the labeling) and remain sterile as long as the package integrity has not been violated. Inspect each package prior to use and do not use the component if any seal or cavity is damaged or breached or if the expiration date has been exceeded. Once opened, the component must be used or discarded.

## **RESTERILIZATION INFORMATION**

Resterilization of this device is not recommended.

## **PATIENT COUNSELING INFORMATION**

Complications and/or failure of total hip prostheses are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Physical activity can result in loosening, wear, and/or fracture of the hip implant. The prospective implant patient must be counseled about the capabilities of the implant and the impact it will have on his or her lifestyle. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the

possibility that the implant or its components may wear out, fail, or need to be replaced. The implant may not last the rest of the patient's life, or any specified length of time. Because prosthetic joints are not as strong, reliable, or durable as a natural, healthy joint, all prosthetic hips may need to be replaced at some point.

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