87-6203-883-99

**R<sub>x</sub> only**

Rev. D September 2017  
© 2017, 2013, 2011 Zimmer, Inc.

## NEXGEN® CR-FLEX AND LPS-FLEX FIXED BEARING KNEE

**zimmer**

Zimmer, Inc.  
1800 West Center Street  
Warsaw, Indiana 46580  
USA  
1-800-348-2759 (US only)  
+1-574-372-4999

Zimmer U.K. Ltd.  
9 Lancaster Place  
South Marston Park  
Swindon, SN3 4FP, UK

### ENGLISH

The CE mark is valid only if it is also printed on the product label.

### NEXGEN® CR-FLEX AND LPS-FLEX FIXED BEARING KNEE

Before using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g., product literature, written surgical technique). Zimmer is not liable for complications that may arise from use of the device in circumstances outside of Zimmer's control including, but not limited to, product selection and deviations from the device's indicated uses or surgical technique.

## **DESCRIPTION**

The NexGen® Flex fixed bearing system consists of two styles of semiconstrained, nonlinked condylar knee prostheses:

- CR-Flex Fixed (CR-Flex femoral, NexGen CR articular surface and a NexGen tibial baseplate)
- LPS-Flex Fixed (LPS-Flex femoral, LPS-Flex articular surface and a NexGen tibial baseplate)

The CR-Flex and LPS-Flex femoral components are designed to have a maximum active flexion of 155 degrees. Femoral components are made from Zimaloy® Cobalt-Chromium-Molybdenum alloy or Tivanium® Ti-6Al-4V alloy (check product label for material composition), articular surface and patella components from ultra-high molecular-weight polyethylene and NexGen tibial baseplates are made from Tivanium Ti-6Al-4V alloy. See the separate package inserts for primary porous and augmentation patella information. A variety of stem extensions, augments, and wedges are available. Articular surfaces are available in multiple thicknesses to facilitate soft tissue tensioning and joint line restoration.

NexGen components are manufactured in a variety of conditions (i.e., porous, precoat, and option [noncoated]). Therefore, some femoral and tibial baseplate components include poly(methyl methacrylate) [PMMA] or hydroxyapatite/tricalcium phosphate (HA/TCP)\*. Please call your Zimmer representative for specific design configuration and surface coating availability.

## **INDICATIONS**

- This device is indicated for patients with severe knee pain and disability due to:
  - Rheumatoid arthritis, osteoarthritis, traumatic arthritis, polyarthritis.
  - Collagen disorders, and/or avascular necrosis of the femoral condyle.
  - Post-traumatic loss of joint configuration, particularly when there is patellofemoral erosion, dysfunction or prior patellectomy.
  - Moderate valgus, varus, or flexion deformities.
  - The salvage of previously failed surgical attempts or for a knee in which satisfactory stability in flexion cannot be obtained at the time of surgery.

- CR-Flex and LPS-Flex porous coated femoral and tibial baseplate components may be used cemented or uncemented (biological fixation). CR-Flex HA/TCP\* porous coated femoral and tibial baseplate components may only be used uncemented. All other femoral, tibial baseplate and all-polyethylene patella components are indicated for cemented use only.
- Specific uses with CR-Flex or LPS-Flex femorals:
  - Provides increased flexion capability for patients who have both the flexibility and desire to increase their flexion range.
  - The CR-Flex femoral, when used with 90-prefix NexGen CR articular surfaces, 00-prefix 10, 12 or 14 mm or 90-prefix 17 or 20 mm Prolong® Highly Crosslinked Polyethylene CR articular surfaces, is designed for use with a functional posterior cruciate ligament and when load bearing range of motion (ROM) is expected to be less than or equal to 155 degrees.
  - The LPS-Flex femoral, when used with the LPS-Flex articular surfaces, is designed for use with both cruciate ligaments excised and when load bearing ROM is expected to be less than or equal to 155 degrees.
- Specific Uses with NexGen LPS or LCCK femorals:
  - The LPS femoral is designed for use with both cruciate ligaments excised and when load bearing ROM is expected to be less than or equal to 130 degrees.
  - The usage of the LCCK femoral is the same as the LPS with the exception that femoral bone loss can be accommodated by femoral stem extensions and augmentation. No varus/valgus constraint or other stability other than that provided by the NexGen LPS is provided when the LCCK femoral component is used with the LPS-Flex articular surfaces.

## **INDIVIDUALIZATION OF TREATMENT**

- Check the appropriate knee implant size matching chart for component matching instructions. Mismatching may result in poor surface contact and could produce pain, decrease wear resistance, produce instability of the implant, or otherwise reduce implant life.
- Use only instruments and provisionals specifically designed for use with these devices to help ensure accurate surgical implantation, soft tissue balancing, and evaluation of knee function.

- Selection of polyethylene components is a matter of physician discretion. Thicker polyethylene components may be needed if the patient is young, heavy, and/or physically active.

## **CONTRAINDICATIONS**

- Contraindications include:
  - Previous history of infection in the affected joint and/or local/systemic infection that may affect the prosthetic joint.
  - Insufficient bone stock on femoral or tibial surfaces.
  - Skeletal immaturity.
  - Neuropathic arthropathy.
  - Osteoporosis or any loss of musculature or neuromuscular disease that compromises the affected limb.
  - A stable, painless arthrodesis in a satisfactory functional position.
  - Severe instability secondary to the absence of collateral ligament integrity.
- Total knee arthroplasty is contraindicated in patients who have rheumatoid arthritis (RA) and an ulcer of the skin or a history of recurrent breakdown of the skin because their risk of postoperative infection is greater. RA patients using steroids may also have increased risk of infection. Late infections in RA patients have been reported 24+ months postoperative.

## **WARNINGS**

- This device is for single patient use only. Do not reuse.
- Do not reinsert an articular surface implant that has been inserted previously. There may be visually nondetectable flaws that could reduce the service life of the implant.
- Prior to closure, the surgical site should be thoroughly cleansed of bone chips, bone cement, and any other debris. Foreign particles at the articular surface may cause excessive wear.
- Do not use without bone cement if bone quality is poor or incapable of providing good initial fixation.
- If uncemented, the patient must be limited in activity and avoid extreme positions that place strain on the knee during the first 12 weeks post-operative.

- Do not use:
  - This product for other than labeled indications (off-label use).
  - Any component if damage is found or caused during setup or insertion.
  - Components from other knee systems (and vice versa) unless expressly labeled for such use. Premature wear or loosening may develop and may require surgical explantation. Exception: NexGen stem extensions are compatible with Insall/Burstein® II Knee tibial and femoral components. Insall/Burstein II stem extensions\* are compatible with NexGen tibial but not NexGen femoral components (they lack a set screw locking groove).
  - LPS-Flex femoral components with any NexGen CR or PS articular surface components. Do not use CR-Flex femoral components with LPS, LPS-Flex, or PS articular surface components. They were not designed to be compatible.
  - LPS or LPS-Flex femoral components with LPS-Flex Prolong articular surfaces unless the part number of the femoral component has a 51 or 52 part number suffix. Use of other suffix LPS and LPS-Flex femoral components may increase the risk of articular surface spine fracture. LPS-Flex Prolong articular surfaces may be used with LCCK femoral components.
  - LPS-Flex Prolong articular surfaces for revision surgery of an LPS or LPS-Flex knee prosthesis unless it can be confirmed that the implanted femoral component has a 51 or 52 part number suffix. If this cannot be confirmed, then replace the femoral component with an LPS or LPS-Flex femoral component with a 51 or 52 part number suffix, or with an LCCK femoral component.
  - LPS-Flex femoral components with LPS articular surface components. They were not designed to be compatible.
  - 4-pegged tibial baseplates with LPS, LPS-Flex or LCCK femoral or 17 or 20 mm 4-pegged tibial baseplates with CR-Flex femoral components. The fixation integrity is unknown.
  - Micro sized patellas with CR-Flex C-H size or LPS-Flex femoral components. Excessive wear may result. For A or B sized CR-Flex femoral components, use only micro sized patellas.
  - 26 mm sized standard patellas with any size CR-Flex or LPS-Flex femoral components unless used in an inset mode. Excessive wear may result.
  - 29 or 32 mm sized standard patellas with size G or H LPS-Flex femoral components unless used in an inset mode. Excessive wear may result.

- CR-Flex femoral components recommended for use with 90-5970 series CR, 00-5952 series 10 – 14 mm Prolong CR or 90-5952 series 17 or 20 mm Prolong CR articular surfaces.
- Use LPS-Flex fixed bearing or NexGen LPS articular surfaces only with LPS femoral components.
- Use LPS-Flex fixed bearing, NexGen LPS or LCCK articular surfaces only with LCCK femoral components.
- Use LPS-Flex Titanium femoral components only with LPS-Flex Prolong articular surfaces.
- Use only NexGen all-polyethylene, porous or augmentation patellas with the CR-Flex, LPS-Flex, LPS or LCCK femoral components. Patellas made for other systems may demonstrate excessive wear when used with NexGen Knee femoral components.
- All CR-Flex and LPS-Flex 17, 20 and 23 mm articular surfaces require a locking screw to fasten the articular surface to the tibial baseplate.
- Screw fixation is required when using uncemented porous 4-pegged or stemmed tibial baseplates.
- Stem extensions are required when a tibial baseplate is being used with an LCCK articular surface.
- The risk of implant failure is higher with inaccurate component alignment or positioning.
- Soft tissues should be balanced and components positioning confirmed to minimize edge loading.
- Fat embolism risk is increased with intramedullary instrumentation and/or cement pressurization. Consider venting the femur or tibia.
- With simultaneous bilateral knee surgery, release leg tourniquets ten minutes apart to lessen any lung insult that may occur.

## **PRECAUTIONS**

- Avoid notching, scratching, or striking the device.
- The potential for deep sepsis can be minimized by using biocontamination controls. Continued surveillance for new or recurrent sources of infection should be continued as long as the device is in place.
- Use of this device may cause patella impingement at high flexion if used with 00-5970 or 00-5966 series NexGen CR articular surfaces.

## **MAGNETIC RESONANCE (MR) SAFETY AND COMPATIBILITY**

In the event of exposure to foreseeable environmental conditions, such as magnetic fields, the user and/or patient should be informed of the following precautions:

- This device has not been evaluated for safety and compatibility in the MR environment.
- This device has not been tested for heating or migration in the MR environment.
- There is a potential for heating and migration in the MR environment.
- There is the potential for metal implants to create MR imaging artifacts in the vicinity of the implant.

## **ADVERSE EFFECTS**

- Loosening or fracture/damage of the prosthetic knee components or surrounding tissues.
- Dislocation and/or joint instability.
- Malalignment of the prosthetic knee components.
- Bone fracture or nerve damage.
- Swelling or infection.
- Leg length discrepancies.
- Poor range of motion.
- Pain.
- Venous thromboembolic disease.
- Inflammation.
- Metal sensitivity.
- Corrosion of metal components (the significance and long-term implications are uncertain and await further clinical evidence and evaluation).
- Wear debris can initiate osteolysis which may result in loosening of the implant.

## STERILITY

These devices are provided sterile (sterilized by Gamma irradiation – indicated by the **STERILE R** symbol on the labeling, sterilized by STERRAD gas plasma sterilization – indicated by the **STERILE GP** symbol on the labeling, or sterilized by ethylene oxide – indicated by the **STERILE EO** symbol on the labeling) and remain sterile as long as the package integrity has not been violated. Inspect each package prior to use and do not use the component if any seal or cavity is damaged or breached or if the expiration date has been exceeded. Once opened, the component must be used, discarded or resterilized.

## RESTERILIZATION INFORMATION

- These sterilization instructions are consistent with ANSI/AAMI/ISO standards and guidelines. They should be used for sterile items that were opened but unused.
- Do not reuse instruments or devices labeled for single use only.
- In the event of inadvertent loss of sterility while preparing for surgery, sterile packaged, all metal implants may be resterilized only once for immediate use. This is subject to the exceptions listed below.
- DO NOT RESTERILIZE:
  - Single use only components that have been contaminated with body fluids or debris or previously implanted.
  - Trabecular Metal™ Technology components.
  - Components with a packaging expiration date that has been exceeded.
  - Components containing UHMWPE or HDPE.
  - Components containing PMMA.
  - Components containing polyaryletherketone (PAEK) or polyetheretherketone (PEEK).
  - Hydroxyapatite/tricalcium phosphate (HA/TCP) coated components.
  - Plastic foam components.

- Do not use the original plastic cavities or lids for resterilization. Single devices should be packaged in a medical grade sterilization pouch or wrap which conforms to the recommended specifications for steam sterilization provided in the table below. Ensure that the pouch or wrap is large enough to contain the devices without stressing the seals or tearing the pouch.
- Rinse porous components to remove lint or debris (using USP purified water).
- Aggressive cleaning with detergents and brushes may damage special features of the implant, such as fiber metal pads or bead coatings. Also, certain detergents may be difficult to rinse off polymer items, especially those made of silicone rubber.
- Items made from titanium and titanium alloys can form oxide layers from steam boiler treatment chemicals or detergent residues. While these oxides are biocompatible, they can obliterate etchings and stampings.
- Modular implant components must be sterilized separately to minimize potential bioburden buildup in the dead space and expansion/contraction stresses.
- Never subject ceramic components to sudden cooling or quenching. Allow them to cool slowly.

### **Recommended Sterilization/Resterilization Specifications**

Follow the sterilizer manufacturer's instructions for loading patterns and selection of sterilization parameters. Drying times vary according to load size and should be increased for larger loads.

## Single, All Metal Implants

### Steam Sterilization

Type	Temperature	Exposure Time	Minimum Dry Time
Gravity Displacement	121°C (250°F)	30 minutes	15 minutes
Gravity Displacement	132°C (270°F)	15 minutes	
UK Pre-vacuum/ Pulsating Vacuum <sup>1</sup>	134°C (273°F)	3 minutes	
Pre-vacuum/ Pulsating Vacuum	132°C (270°F)	4 minutes	

<sup>1</sup> This cycle is not for use in the United States.

Please contact Zimmer at the following number if you have additional questions. In the USA, call 1-800-348-2759. For calls outside the USA, call the local international access code +1-574-267-6131.

### **PATIENT COUNSELING INFORMATION**

Complications and/or failure of total knee prostheses are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients that fail to follow through with the required rehabilitation program. Not all patients with CR-Flex or LPS-Flex components will achieve 155 degrees of flexion. Excessive physical activity and injury can result in loosening, wear, and/or fracture of the knee implant. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail, or need to be replaced. The implant may not, and is not guaranteed to, last the rest of the patient's life. Because prosthetic joints are not as strong, reliable, or durable as natural, healthy joints, all prosthetic knees may need to be replaced at some point.

\*not licensed in Canada

† Trademark of The Hospital for Special Surgery

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