BIOMET

- UK** POLYETHYLENE COMPONENTS
ATTENTION OPERATING SURGEON
- CZ** POLYETHYLENOVÉ KOMPONENTY
UPOZORNĚNÍ PRO OPERUJÍCÍHO CHIRURGA
- DE** POLYETHYLENKOMPONENTEN
INFORMATION FÜR DEN OPERIERENDEN CHIRURGEN
- DK** POLYETHYLENKOMPONENTER
BEMÆRKNING TIL DEN OPERERENDE KIRURG
- EL** ΕΞΑΡΤΗΜΑΤΑ ΠΟΛYΑΙΘYΛΕΝΙΟΥ
ΥΠΟΧΗ ΤΟΥ ΧΕΙΡΟΥΡΓΟΥ
- ES** COMPONENTES DE POLIETILENO
A LA ATENCIÓN DEL CIRUJANO
- FI** POLYETEENIOSAT
LEIKKAVALLE KIRURGILLE
- FR** PIÈCES EN POLYÉTHYLÈNE
À L'ATTENTION DU CHIRURGIEN
- HU** POLIETILÉN RÉSZEK
AZ OPERÁLÓ SEBÉSZ FIGYELMÉBE
- IT** COMPONENTI IN POLIETILENE
ALL'ATTENZIONE DEL CHIRURGO
- KO** 폴리에틸렌 구성 요소
시술 시 주의 요망
- NL** POLYETHYLEEN COMPONENTEN
INFORMATIE VOOR DE BEHANDELEND CHIRURG
- NO** KOMPONENTER AV POLYETYLEN
INFORMASJON FOR KIRURGEN
- PL** KOMPONENTY POLIETYLENOWE
UWAGI DLA OPERUJĄCEGO CHIRURGA
- PT** COMPONENTES EM POLIETILENO
À ATENÇÃO DO CIRURGIÃO
- RU** КОМПОНЕНТЫ ИЗ ПОЛИЭТИЛЕНА
ВНИМАНИЮ ОПЕРИРУЮЩЕГО ХИРУРГА
- SE** POLYETYLENKOMPONENTER
ATT OBSERVERA FÖR OPERERANDE ORTOPED
- TR** POLIETİLEN BİLEŞENLER
CERRAHİN DIKKATİNE
- ZH** 聚乙烯组件
主刀外科医生需知

BIOMET UK LTD.
WATERTON INDUSTRIAL ESTATE
BRIDGEND CF31 3XA, UK
TEL: +44 (0) 1656 655221
FAX: +44 (0) 1656 645454

REF: 5401000258
REVISION: E
DATE: FEB 2017
Change Control No: 229338

For additional information please see www.biomet.com

UK

POLYETHYLENE COMPONENTS

ATTENTION OPERATING SURGEON

DESCRIPTION

Acetabular liners and cemented cups are manufactured from ultra high molecular weight polyethylene and comprise various designs and sizes. The acetabular components are utilised in conjunction with other hip prosthesis components to form part of a total joint replacement system. Total hip joint prostheses may include: femoral stems, femoral heads, acetabular shells, acetabular liners and cups. Biomet has evaluated the compatibility of Biomet devices with implants and components from Zimmer orthopaedic companies. Only authorized combinations should be used. To determine whether these devices have been authorized for use in a proposed combination with Zimmer products, please contact your sales representative and reference the following website: <http://www.productcompatibility.zimmer.com>

MATERIALS

| | |
|----------------|--|
| ArCom/UHMWPE | Ultra high molecular weight polyethylene |
| ArComXL | Highly cross linked ultra high molecular weight polyethylene |
| E1 | Highly cross linked ultra high molecular weight polyethylene with α -tocopherol |
| Cement Spacers | Polymethylmethacrylate |
| Marker Wire | Titanium-aluminium-vanadium alloy |

INDICATIONS

The indications for the use of acetabular liners and cemented cups include:

1. Non-inflammatory degenerative joint diseases including osteoarthritis, avascular necrosis and post-traumatic arthritis.
2. Rheumatoid arthritis.
3. Correction of functional deformity.
4. Revision of failed joint reconstruction or treatment.
5. Treatment of femoral neck and trochanteric fractures of the proximal femur with femoral head involvement and which are unmanageable using other techniques.

CONTRAINDICATIONS

Absolute contraindications: infection, sepsis and osteomyelitis

Relative contraindications:

1. Uncooperative patient or patient with neurologic disorders who are incapable of following instructions
2. Osteoporosis
3. Metabolic disorder which may impair bone formation
4. Osteomalacia
5. Local and distant foci of infection
6. Rapid joint destruction, marked bone loss or bone resorption apparent on roentgenogram
7. Vascular insufficiency, muscular atrophy, or neuromuscular disease
8. Skeletal immaturity
9. Morbid obesity
10. Foreign body sensitivity. Where suspected or unknown, material sensitivity tests are to be made prior to implantation
11. Any condition that may interfere with the survival of the implants such as Paget's disease, Charcot's disease, sickle cell anaemia or traits, lower extremity muscular atrophy or neuromuscular disease

PATIENT SELECTION FACTORS

Patient selection factors which should be considered include:

1. Primary need to obtain pain relief.
2. Primary need to improve function.
3. Ability of the patient to follow instructions, including control of weight and activity level, and a willingness to undergo rehabilitation.
4. A good nutritional state of the patient.
5. Whether the patient has reached full skeletal maturity.

WARNINGS AND PRECAUTIONS

1. Improper selection, placement, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may lead to subsequent reduction in the service life of the prosthetic components.
2. Malalignment of components or inaccurate implantation can lead to excessive wear and/or failure of the implant or procedure.
3. Complete preclosure cleaning and removal of bone cement debris, metallic debris and other surgical debris at the implant site is critical to minimise wear of the implant articular surfaces.
4. Improper preoperative or intra-operative implant handling or damage (scratches, dents, etc.) can lead to crevice corrosion, fretting, fatigue fracture, and/or excessive wear.
5. Do not modify implants unless specified in the operative technique.
6. The surgeon is to be thoroughly familiar with the implants, instruments and surgical procedure prior to performing surgery.
7. Biomet joint replacement prostheses provide the surgeon with a means of reducing pain and restoring function for many patients. While these devices are generally successful in attaining these goals, they cannot be expected to withstand the same activity levels and loads of normal healthy bone and joint tissue.

- Specialised instruments are designed for specific Biomet hip systems to assist in the accurate implantation of the implants. The use of other manufacturer's instruments can result in inaccurate fit and sizing.
- Whilst rare, intraoperative fracture or breakage of instruments has been reported. Instruments that have experienced excessive use or excessive force are susceptible to fracture. Biomet recommends that instruments are examined for wear or damage prior to surgery.
- Do not reuse implants. While an implant may appear undamaged, previous stress may have created imperfections that would reduce the service life of the implant. Do not treat patients with implants that have been, even momentarily, placed in a different patient. Furthermore, re-using an implant could cause patient contamination.
- Use clean surgical gloves when handling implants.
- Do not use any component that is nicked, scratched, chipped or otherwise altered.
- Do not use prosthetic implants from other systems with Biomet components due to the probability of incompatible sizing and bearing surfaces, which may lead to premature wear, malalignment and failure. For specific Zimmer/Biomet product compatibility, reference the following website: www.productcompatibility.zimmer.com
- All modular components must be accurately seated to minimize the risk of disassociation.
- Prior to seating the liner into the shell component, all surgical debris (tissue fragments, etc.) must be removed from the interior of the shell component, as debris may inhibit the locking mechanism from engaging and securing the liner into the shell component.
- In any instance where a liner engages the locking ring and the liner is subsequently removed or replaced, the locking ring should be replaced with a new ring.
- Bone stock of adequate quality must be present and appraised at the time of surgery.
- Patients should be warned of the impact of excessive loading that can result if the patient is involved in an occupation that includes substantial walking, running, lifting, or excessive muscle loading due to weight that place extreme demands on the hip and can result in device failure or dislocation.
- Patient smoking may result in delayed healing, non-healing and/or compromised stability in or around the placement site.
- Deviations from the recommended technique could result in early loosening/failure of the device, or other adverse events as outlined in the following section.

Specific warnings for cemented fixation:

- Care must be taken to ensure complete support of all parts of the device embedded in bone cement to reduce the risk of stress concentrations which may lead to failure of the implant.
- Complete removal of all bone cement debris is critical to minimise or reduce the risk of wear of the articulating surfaces.
- Implant fracture due to cement failure has been reported.
- Please refer to IFU supplied with cement for further information.

Please refer to the operative technique guide available from your Biomet distributor for suggested operative technique.

MRI COMPATIBILITY

The risks associated with these devices have been evaluated for safety and compatibility in the MR environment. The risks are known to include heating, migration, and image artefacts at or near the implant site. Scanning a patient who has this device may result in patient injury.

PATIENT WARNINGS

Postoperative care is important. The patient must be warned of the limitations of the reconstruction and the need for protection of the implants from full weight bearing or load bearing until adequate fixation and healing has occurred. The patient is to be cautioned to govern activities, protecting the joint replacement from unreasonable stress conditions. Excessive activity, failure to control body weight and trauma affecting the joint replacement have been associated with premature failure of the reconstruction by loosening, fracture and/or wear of the implants. Loosening of the components can result in increased production of wear particles, as well as accelerating damage to bone, making successful revision surgery more difficult. The patient is to be advised of the importance of postoperative follow-up examination and that excessive activity or trauma can lead to failure. The patient is to be warned of the surgical risks and made aware of possible adverse effects in advance of surgery. The patient must be warned that the device does not replace normal healthy bone, and that the implant can break or be damaged as a result of excessive load bearing or trauma. The patient must be warned to inform any other medical practitioner who may treat him in the future of the presence of the implant.

POSSIBLE ADVERSE EFFECTS

- Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions involving various sizes of macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to or during the healing process. Particulate wear debris and discoloration from metallic and polyethylene components of joint implants may be present in adjacent tissue or fluid. It has been reported that wear debris may initiate a cellular response resulting in osteolysis or osteolysis may be a result of loosening of the implant. An association had been reported between articulating surfaces of: i) CoCrMo alloy on polyethylene, and ii) Titanium alloy on polyethylene in hip replacements and increased genotoxicity. However the clinical relevance of the data is not known and no definite conclusions were drawn as to which metal ions or interactions between metal ions or particulate metals might be responsible for the observed data. It was also noted that an association does not necessarily mean a causal relationship, and that any potentially increased risk associated with metal ions needs to be balanced against the benefits resulting from hip replacement.
- Early or late postoperative infection and allergic reaction.
- Intraoperative bone perforation or fracture may occur, particularly in the presence of poor bone stock caused by osteoporosis, bone defects from previous surgery, bone resorption, or while inserting the device.
- Loosening or migration of the implants can occur due to loss of fixation, trauma, malalignment, bone resorption, and/or excessive activity.
- Periarticular calcification or ossification with or without impediment of joint mobility.
- Inadequate range of motion due to improper selection or positioning of components.
- Leg length discrepancy.
- Dislocation and subluxation due to inadequate fixation and improper positioning. Muscle and fibrous tissue laxity can also contribute to these conditions.
- Fatigue fracture of component can occur as a result of loss of fixation, strenuous activity, malalignment, trauma, non-union, and/or excessive weight.
- Fretting and crevice corrosion can occur at interfaces between components.
- Wear and/or deformation of articulating surfaces.

12. Trochanteric avulsion or non-union as a result of excess muscular tension, early weight bearing, or inadequate reattachment.
13. Problems of the knee or ankle of the affected limb or contralateral limb aggravated by leg length discrepancy, too much femoral medialisation or muscle deficiencies.
14. Intraoperative or postoperative bone fracture and/or postoperative pain.
15. Peripheral neuropathies have been reported following total joint replacement surgery. Subclinical nerve damage occurs more frequently, possibly as a result of surgical trauma.
16. Cardiovascular reactions: venous thrombosis, pulmonary embolism or myocardial infarction.
17. Haematoma.
18. Delayed wound healing.
19. Altered gait / limp.
20. Bursitis.
21. Gastrointestinal or genitourinary complications.
22. Pulmonary embolism.
23. Death.
24. Effusion.

STERILITY

Unless otherwise specified Biomet joint replacement components are sterilised by gamma irradiation (25-40 kGy). This is denoted on the label by the symbols **STERILE**.

ArcomXL Highly Crosslinked UHMWPE acetabular components are sterilised by exposure to a Gas Plasma sterilisation process. This is denoted on the label by the symbol **STERILE**.

The packages should be inspected for damage or flaws in the sterile barrier before opening. Product should be assumed to be non-sterile if the packaging is damaged. Do not use if packaging is damaged.

CE Mark on the package insert (IFU) is not valid unless there is a CE Mark on the product (description) label.

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