## ENGLISH

The CE mark is valid only if it is also printed on the product label.

### **Bone Screws for Acetabular Cup**

#### **Important information for the operating surgeon**

Before using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g., product literature, written surgical technique). Zimmer is not liable for complications arising from the use of the device outside of its indicated uses, surgical technique or judgment, product selection, and similar matters outside the control of Zimmer.

#### **Product Compatibility**

Due to the acquisition of pre-existing product lines, Zimmer has initiated a testing program to evaluate the compatibility of these devices with implants and components made or distributed by all Zimmer orthopaedic companies, which include Zimmer GmbH (previously Centerpulse Orthopedics Ltd.), Zimmer, Inc., Zimmer Trabecular Metal Technology, Inc. (previously Implex Corp.), Zimmer U.K. Ltd., and Zimmer Austin, Inc. (previously Centerpulse Orthopedics, Inc.). Only authorized combinations should be used. To determine whether these devices have been authorized for use in a proposed combination, contact your Zimmer sales representative or visit the Zimmer website: [www.productcompatibility.zimmer.com](http://www.productcompatibility.zimmer.com). A printout of the website information can also be obtained by calling Zimmer, Inc. Customer Service, 1-800-348-2759 (U.S.) or the local international access code +1-574-372-4999 (outside the U.S.). Former Centerpulse and Implex products that are now packaged in Zimmer boxes, and for which compatibility could be an issue, have been labeled "former Centerpulse" and "former Implex" to provide clarification for the user.

#### **This Physicians Insert is valid for the following components identified below:**

- **Cedior™ Universal screw (*Protasul*®-WF64, Ti-6Al-4V)**
- **Countersunk Cancellous Bone Screw (*Protasul*-100, Ti-6Al-7Nb)**
- **Countersunk Cancellous Bone Screw (*Protasul*-WF64, Ti-6Al-4V)**

## **DESCRIPTION**

The Bone Screw is used in conjunction with an acetabular cup component for replacement of the acetabulum in total hip arthroplasty. The Bone Screws are available in different sizes and different lengths. Zimmer implants for internal fixation are available in the following titanium versions, depending on the product.

Titanium Alloy *Protasul-64* (Ti-6Al-4V) or Titanium Alloy *Protasul-100* (Ti-6Al-7Nb) are used.

For detailed information refer to the corresponding Instructions for Use of the implants the screws should be used with.

## **INDICATIONS/INTENDED USE**

Bone screws for acetabular cup are used to enable adequate fixation at specified acetabular implants in the course of total hip arthroplasty, thereby helping to reduce pain and increase hip mobility through cementless fixation in the acetabulum and ilium of patients with adequate bone stock to support the component.

The indications for the respective Bone Screws for acetabular cup are determined by the acetabular implant used in conjunction with them. For Cages and Rings, the supplementary fixation of the implant with screws is mandatory, whereas for the other shells the use of screws for additional fixation is optional.

## **CONTRAINDICATIONS**

- All concomitant diseases that may impair the fixation of the implant and/or the success of the intervention.
- Lack of bone substance or poor bone quality impairing stable fixation of the implant.
- Allergy to the implanted material.
- Acute or chronic, local or systemic infections.
- Further product dependent contraindications are mentioned in the corresponding Instructions for Use of the implants the screws should be used with.

## **WARNINGS**

- Implants are for single use only. Do not reuse.
- Do not use any component if damage is found or caused during setup or insertion.
- Implants and implant parts must only be combined with components belonging to the same system. No liability is accepted for products of third parties.

- Use only instruments specifically designed for use with these devices to help ensure accurate surgical implantation.
- Complications or failure of any total hip prosthesis are more likely to occur in heavy patients.
- The load-bearing capacity of the implant can be compromised by notching, scratching, or striking.
- Improper selection, placement, positioning, and fixation of the implant components may result in unusual stress conditions reducing the service life of the implants.
- Do not use these products for other than labeled indications (misuse).

### **PRECAUTIONS**

- Continued surveillance for new or recurrent infection should be continued as long as the device is in place.
- The risk of adverse effects such as secondary infections, allergies, material fatigue (breakage), implant failure and/or impaired blood circulation, increases with duration of time *in situ*, non-union, patient weight and activity level.
- Implants must not be machined or altered in any way, unless this is expressly envisaged in the design and in the surgical technique.

### **ADVERSE EFFECTS**

- |   |   |
|---|---|
| • Bone fractures  | • Hypersensitivity to implant material  |
| • Vascular complications                                      | • Inflammatory reactions and osteolysis |
| • Corrosion of metal implants                                 | • Nerve damage                          |
| • Damage of surrounding tissues                               | • Pain                                  |
| • Disassociation of components                                | • Perforation of the acetabulum         |
| • Displacement and loosening of implant or implant components | • Wear                                  |
| • Fracture of the implant                                     | • Wound infections                      |

### **STERILIZATION INSTRUCTIONS**

- These sterilization instructions are consistent with ANSI/AAMI/ISO standards and guidelines. They should be used for items supplied non-sterile, for reprocessing reusable devices, or for sterile items that were opened but unused.

- See Zimmer booklet 97-5000-170-00 (available from your distributor or at [www.zimmer.com](http://www.zimmer.com)) for detailed guidelines.
- Do not reuse instruments or devices labeled for single use only. All metal implants without an articulating surface may be resterilized multiple times. This is subject to the exceptions listed below.
- DO NOT RESTERILIZE:
  - Single use only components that have been contaminated with body fluids or debris or previously implanted.
  - Components with a packaging expiration date that has been exceeded.
  - Components containing UHMWPE or HDPE.
- Do not use the original plastic cavities or lids for resterilization. Single devices should be packaged in a medical grade sterilization pouch or wrap which conforms to the recommended specifications for steam sterilization provided in the table below. Ensure that the pouch or wrap is large enough to contain the devices without stressing the seals or tearing.
- Sets of devices may be loaded into general or dedicated trays. Use a medical grade sterilization pouch or wrap (double wrap method – ANSI/AAMI ST79). Sterilization should be performed immediately prior to use.
- Items made from titanium and titanium alloys can form oxide layers from steam boiler treatment chemicals or detergent residues. While these oxides are biocompatible, they can obliterate etchings and stampings.

### **Recommended Sterilization/Resterilization Specifications**

Follow the sterilizer manufacturer's instructions for loading patterns and selection of sterilization parameters. Drying times vary according to load size and should be increased for larger loads.

## Single, All Metal implants

### **Steam Sterilization**

Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Gravity Displacement	121°C (250°F)	30 minutes	15 minutes
Gravity Displacement	132°C (270°F)	15 minutes	
UK Pre-vacuum/ Pulsating Vacuum <sup>1</sup>	134°C (273°F)	3 minutes	
Pre-vacuum/Pulsating Vacuum	132°C (270°F)	4 minutes	

<sup>1</sup> This cycle is not for use in the United States.

## All Metal Implants In Sterilization Cases and Trays

### **Steam Sterilization**

Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
UK Pre-vacuum/ Pulsating Vacuum <sup>1</sup>	134°C (273°F)	3 minutes	30 minutes
Pre-vacuum/Pulsating Vacuum	132°C (270°F)	4 minutes	
Pre-vacuum/Pulsating Vacuum <sup>1,2</sup>	134°C (273°F)	18 minutes	

<sup>1</sup> This cycle is not for use in the United States.

<sup>2</sup> This cycle is not to be used for the inactivation of prions.

Please contact Zimmer at the following number if you have additional questions. In the USA, call 1-800-348-2759. For calls outside the USA, call the local international access code +1-574-267-6131.

### **STORAGE AND HANDLING**

- Implants, implant parts and instruments that can no longer be used may be returned to the manufacturer for proper disposal free of charge.

### **PATIENT COUNSELING INFORMATION**

Complications and/or failure of prosthetic implants are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Physical activity or trauma can result in loosening, wear, and/or fracture of the implant. The patient must be counseled about the capabilities of the implant and the impact it will have on his or her lifestyle. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant may not last the rest of the patient's life, or any particular length of time. Because prosthetic implants are not as strong, reliable, or durable as natural, healthy tissues/bones, all such devices may need to be replaced at some point.

Caution: Federal law restricts this device to sale by or on the order of a physician

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