87-6203-826-99	<b>R<sub>x</sub> only</b>		Rev. B January 2017 © 2017, 2012, Zimmer, Inc.
<b>VIVACIT-E<sup>®</sup> VITAMIN E HIGHLY CROSSLINKED POLYETHYLENE LINERS</b>			
 <b>zimmer</b>	 Zimmer, Inc. 1800 West Center Street Warsaw, Indiana 46580 USA 1-800-348-2759 (US only) +1-574-372-4999	 Zimmer U.K. Ltd. 9 Lancaster Place South Marston Park Swindon, SN3 4FP, UK	 0086
			

## ENGLISH

The CE mark is valid only if it is also printed on the product label.

## VIVACIT-E<sup>®</sup> VITAMIN E HIGHLY CROSSLINKED POLYETHYLENE LINERS

Before using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g., product literature, written surgical technique). Zimmer is not liable for complications that may arise from the use of the device in circumstances outside of Zimmer's control including, but not limited to, product selection and deviations from the device's indicated uses or surgical technique.

Due to the acquisition of pre-existing product lines, Zimmer has initiated a testing program to evaluate the compatibility of these devices with implants and components made or distributed by all Zimmer orthopaedic companies, which include Zimmer GmbH (previously Centerpulse Orthopedics Ltd.), Zimmer, Inc., Zimmer Trabecular Metal Technology, Inc. (previously Implex Corp.), Zimmer U.K. Ltd., and Zimmer Austin, Inc. (previously Centerpulse Orthopedics, Inc.) and Biomet Orthopedics companies. Only authorized combinations should be used. To determine whether these devices have been authorized for use in a proposed combination, please contact your Zimmer sales representative or visit the Zimmer website: [www.productcompatibility.zimmer.com](http://www.productcompatibility.zimmer.com). A printout of the website information can also be obtained by calling Zimmer, Inc. Customer Service, 1-800-348-2759 (U.S.) or the local international access code +1-574-372-4999 (outside the U.S.). Former Centerpulse and Implex products that are now packaged in Zimmer boxes, and for which compatibility could be an issue, have been labeled “former Centerpulse” and “former Implex” to provide clarification for the user.

## **DESCRIPTION**

The Vivacit-E Vitamin E Highly Crosslinked Polyethylene liners are a series of acetabular liners manufactured from Vitamin E ( $\alpha$ -tocopherol) blended, highly crosslinked ultra-high molecular weight polyethylene (HXPE) and are designed for use with the Continuum<sup>®</sup>, Trilogy<sup>®</sup> IT, and Allofit<sup>®</sup> IT Alloclassic<sup>®\*\*\*</sup> acetabular systems.

\*\*Allofit IT Alloclassic Acetabular System is not available in the United States.

The liners are offered in a number of articulating diameters and configurations. Each liner is identified with a two letter size code (e.g. JJ) that corresponds to the mating acetabular shell. Liners are secured to the respective shell by engaging a protruding polyethylene rim with an integrated locking groove in the shell. The Vivacit-E Vitamin E Highly Crosslinked Polyethylene liners are to be used with compatible femoral heads with corresponding articulating diameters.

## **INDICATIONS / INTENDED USE**

- The system is indicated for primary or revision surgery in skeletally mature individuals for rehabilitating hips damaged as a result of noninflammatory degenerative joint disease (NIDJD) or its composite diagnoses of osteoarthritis, avascular necrosis, *protrusio acetabuli*, traumatic arthritis, slipped capital epiphysis, fused hip, fracture of the pelvis, and diastrophic variant.
- The system is intended for use either with or without bone cement in total hip arthroplasty.

## **CONTRAINDICATIONS**

This device is contraindicated for the following:

- Osteoradionecrosis.
- Neuromuscular compromise, vascular deficiency or other conditions in the affected limb that may lead to inadequate skeletal fixation.
- Systemic or local infection.

## **WARNINGS**

- This device is for single patient use only. Do not reuse. Reuse may adversely affect performance of this device and compromise patient safety.
- Do not use any component if damage is found or caused during setup or insertion. Imperfections may reduce the service life of the implant.
- Do not use elevated rim or large inner diameter (ID) liners as a solution for a malpositioned shell. Unsupported loading of polyethylene liners fatigues the polyethylene and may lead to fracturing.
- Improper selection, placement, positioning, and fixation of the implant components may result in unusual stress conditions reducing the service life of the prosthetic implants.
- Do not assemble mating components without ensuring that the surfaces are free of blood and debris. Failure to ensure that mating surfaces are clean and dry could result in inadequate seating of one component upon the other and the subsequent disassembly of the mated components or fracture of the implant.

- Repeated assembly and disassembly of modular components could compromise the critical locking action of the device. Use the provisional components during trial reductions. Change the components only when clinically necessary. If a Vivacit-E polyethylene liner has been engaged and removed for any reason, a new liner must be used.

## **PRECAUTIONS**

- Use only instruments and provisionals designed for use with these devices to help ensure accurate surgical implantation, soft tissue balancing, and evaluation of joint function.
- Do not use the femoral head or acetabular liner implant components to perform a trial reduction of the hip. Mate only provisional heads with provisional acetabular liners. Contact of provisional heads and liners with implants can result in abrasion. The provisional liners may be used with the acetabular shell implant.
- Appropriate placement of the components helps to avoid unintended loading of the components and may improve service life of the implants. The surgical technique is available through a Zimmer representative to provide additional detailed information regarding the recommended procedure.
- It is important that components be at room temperature before assembly and implantation. Components stored in a warm environment may be difficult to assemble.
- The larger ID liners (36 and 40mm) are intended to be used in patients that may require additional stability as determined by the physician. These liners are only available in the larger OD sizes and should only be used in patients that have an adequately sized acetabulum.
- Surveillance for new or recurrent sources of infection should be continued as long as the device is in place.
- Implants should only be rinsed with sterile solution such as USP purified water or Ringer's solution.

## ADVERSE EFFECTS

The following complications have been reported:

- Peripheral neuropathies
- Wound infections
- Perforation of the acetabulum or femur
- Heterotopic bone formation
- Wear
- Inflammatory reactions and osteolysis
- Implant fracture
- Pain
- Dislocation or subluxation
- Trochanteric problems
- Nerve damage
- Pelvic or femoral fractures
- Early or late loosening of components of total hip arthroplasty
- Disassociation of modular components
- Cardiovascular and circulatory disorders

Although there is no conclusive evidence of the relationship between orthopaedic implants and malignant tumours, any condition that causes chronic damage to tissues may be oncogenic.

## MAGNETIC RESONANCE (MR) SAFETY AND COMPATIBILITY

Non-clinical testing has demonstrated the Vivacit-E Vitamin E Highly Crosslinked Polyethylene Liners when used with the Continuum, Trilogy IT, and Allofit IT Alloclassic \*\* Acetabular System are MR Conditional and

carries the symbol  on the labeling. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5 and 3.0 Tesla (T).
- Maximum spatial gradient field of 2100 Gauss/cm when used with Stainless Steel hip stem and 2500 Gauss/cm when used with Cobalt-Chromium alloy or Titanium alloy hip stem.
- For landmarks below the umbilicus, maximum whole body averaged specific absorption rate (SAR) of 1 W/kg.
- For landmarks above the umbilicus, maximum whole body averaged SAR of 2 W/kg.
- Normal mode operation of the MR system.
- Legs of patient not touching.

In non-clinical 128 MHz (3.0T) testing, the Continuum, Trilogy IT, and Allofit IT Alloclassic \*\* Acetabular System did not produce any magnetically induced displacement force or torque that would result in migration of the device in the spatial gradient fields specified above.

In non-clinical 64 MHz (1.5T) testing, the Continuum, Trilogy IT, and Allofit IT Alloclassic \*\* Acetabular System produced a temperature rise of less than 3°C at a local background SAR of 1 W/kg, as assessed by reference implant method for 15 minutes of MR scanning in a model number 46-258170G General Electric Signa RF coil. Computer simulation of the electromagnetic fields in the body during 1.5T MRI yielded an estimated *in-vivo* temperature rise of less than 1.6°C for a whole body average SAR of 1 W/kg for a patient landmark in the upper thigh.

In non-clinical 128 MHz (3.0T) testing, the Continuum, Trilogy IT, and Allofit IT Alloclassic \*\* Acetabular System produced a temperature rise of less than 2°C at a local background SAR of 1 W/kg, as assessed by reference implant method for 15 minutes of MR scanning in a General Electric Signa HDx 3T MR system with software version 15/LX/MR release 15.0.M4.0910.a. Computer simulation of the electromagnetic fields in the body during 3.0T MRI yielded an estimated *in-vivo* temperature rise of less than 1.9°C for a whole body average SAR of 1 W/kg for a patient landmark in the upper thigh.

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant. Image artifact was evaluated in a phantom model. The greatest artifact occurred for the gradient echo sequence and extended about 60mm from the device (CoCrMo alloy components) and 40mm from the device (Titanium alloy components). Maximal distortion for the spin echo sequences ranged from 20-30mm. The tests were performed in a General Electric Signa HDx 3T MR system with software version 15/LX/MR release 15.0.M4.0910.a.

## **STERILITY**

Vivacit-E Vitamin E Highly Crosslinked Polyethylene implants are sterilized with Ethylene Oxide Gas and carry the symbol **STERILE EO** on the labeling. These devices remain sterile as long as the package integrity has not been violated. Inspect each package prior to use and do not use the component if any seal or cavity is damaged or breached or if the expiration date has been exceeded. Once opened, the component must be used or discarded.

## **STERILIZATION INSTRUCTIONS**

Resterilization of this device is specifically not recommended.

## **PATIENT COUNSELING INFORMATION**

Complications and/or failure of total hip prostheses are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Physical activity can result in loosening, wear, and/or fracture of the hip implant. The prospective implant patient must be counseled about the capabilities of the implant and the impact it will have on his or her lifestyle. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail, or need to be replaced. The implant may not, and is not guaranteed to, last the rest of the patient's life, or any specified length of time. Because prosthetic joints are not as strong, reliable, or durable as a natural, healthy joint, all prosthetic hips may need to be replaced at some point.

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