
BIOMET®

## Calcibon® Inject

EN	Instructions for Use
DE	Gebrauchsinformation
FR	Mode d'emploi
ES	Modo de Empleo
IT	Istruzioni per l'uso
PT	Instruções para uso
NL	Gebruikersinformatie
SV	Användarinformation
NO	Bruksanvisning
DA	Brugsanvisning
FI	Käyttöohje
CS	Návod k použití
PL	Instrukcja użytkowania
TR	Kullanım Kılavuzu
EL	ΠΛΗΡΟΦΟΡΙΕΣ ΓΙΑ ΤΟΝ ΧΡΗΣΤΗ

## Calcibon® Inject

### Description

Calcibon Inject is a safe and fully bio-compatible material which is designed to act as an osteoconductive medium to support the ingrowth and fusion of adjacent viable bone when placed in an osseous environment.

Calcibon Inject will fully resorb and be replaced by natural bone due to its structure which is essentially identical to the mineral phase of human bone. The hydroxyapatite (HA) Paste component will resorb readily, whilst the hydroxyapatite/ tricalcium phosphate granules will resorb more slowly providing a longer term scaffold support.

Calcibon Inject can in many cases spare the patient the trauma of autograft harvesting and provides a sterile fully synthetic alternative to human or animal cadaver bone, and is free of the risk of disease transmission.

### Composition

Calcibon Inject is a combination device consisting of 83.5% HA Paste and 16.5% bone graft granules (0.5-1.0 mm). HA Paste is composed of nominally 38% hydroxyapatite in water, and porous granules are composed of nominally 60% hydroxyapatite and 40% tricalcium phosphate. Hydroxyapatite and tricalcium phosphate have been subjected to extensive and successful clinical studies for many years.

### Indications

Calcibon Inject is a synthetic osteoconductive bone paste that is intended

to be used for filling and supporting the ingrowth of adjacent viable bone in bone defects that are not intrinsic to the stability of the bone structure.

These defects may be located in the extremities, pelvis or dental, periodontal or mandibular areas.

Calcibon Inject may be used for filling of bone defects in the case of fractures, resection of benign tumours or cysts, and for filling autograft harvest sites.

Calcibon Inject may also be used in defects of the alveolar ridge or wall inclusive extraction sockets, and for filling of periodontal bone defects.

Calcibon Inject may be gently pressed into the defect by hand and moulded to more accurately fit the defect if and as required, and may also be mixed with bone marrow aspirate prior to implantation. Calcibon Inject will act as a temporary scaffold and will remain as a soft paste and is not intended to provide structural support during the healing process. The implant is biocompatible, radio-opaque and will resorb in a controlled way.

### Contraindications

Calcibon Inject is not designed or sold for any other use except as indicated. Do not use Calcibon Inject in the presence of any contraindication. Calcibon Inject is contraindicated where the device is intended to provide structural support in the skeletal system and **MUST NOT BE USED WHERE THE IMPLANTATION SITE IS UNSTABLE AND NOT RIGIDLY FIXATED, OR BE USED TO GAIN SCREW FIXATION.**

**Other contraindications include:**

- Existing acute or chronic infections, particularly at the implantation site
- Severe vascular or neurological disease
- Poorly vascularised implantation site
- Uncontrolled diabetes
- Severe degenerative disease
- Connective tissue collagen disease
- Hypercalcemia, abnormal calcium metabolism
- Inflammatory bone disease
- Malignant tumours
- Severely impaired renal function
- Open epiphyseal plates in paediatric patients

**Precautions**

Calcibon Inject is intended for professional use only, and is only intended for use by surgeons familiar with and skilled in the techniques of bone repair and replacement.

Calcibon Inject is not intended for load bearing applications. It is important to ensure that the area around the implantation site be secured mechanically with rigid fixation or cage to provide structural support and maintain the implant in a static, load free environment. Calcibon Inject must not be used to gain screw fixation.

It is important to maximise the contact between existing bone and the implant to ensure proper bone regeneration.

Calcibon Inject should only be implanted into fresh vital bone to ensure rapid revascularisation.

As with any material, care should be taken to avoid the occurrence of emboli, therefore the highly pressurised

application of Calcibon Inject into a tightly confined defect with venous or arterial access should be avoided.

The effect of Calcibon Inject on patients with the following conditions is not known:

- Long term infection
- Metabolic bone diseases
- Radiation bone therapy
- Cardiovascular diseases

The effects of Calcibon Inject in paediatric patients and the effects of mixing with other substances are unknown.

However no special precautions have been identified at the time of issue.

**Possible Complications**

A successful result may not be achieved in every case.

A secondary operation to remove or replace an implant may be necessary due to surgical error, specific medical conditions or device failure.

**Adverse Events**

Possible adverse effects may include but are not limited to:

- Displacement of the Calcibon Inject due to a load being applied
- Bone deformity or non-viable bone at the implantation site
- Post operative complications including infection, haematoma, oedema, swelling, and fluid accumulation, tissue thinning, and other complications that are possible with surgery.
- Allergic reaction to the product

**Warnings**

The product is gamma sterilized and supplied in sealed double pouches. DO NOT USE IF OPENED, PUNCTURED

OR PRODUCT AND / OR APPLICATOR APPEAR DAMAGED OR PASTE HAS DRIED. Read the expiration date before use and DO NOT USE BEYOND THE EXPIRATION DATE.

It is recommended not to open the inner pouch until a few minutes prior to implantation. Calcibon Inject is for SINGLE USE ONLY. Do not attempt to re-sterilise or re-use due to the risk of cross infection. Calcibon Inject should only be used in procedures where the implant can be adequately contained. Calcibon Inject is opaque to x-rays and may hide areas under or above the implant on a radiograph.

### Application

#### Step 1:

Open both the outer and inner pouches. Remove the cap from the applicator. Calcibon Inject is supplied as a paste in the applicator, ready for use.

#### Step 2:

Implant the quantity of Calcibon Inject as required by the defect. The maximum dose of Calcibon Inject is 20g per defect. The defect should be completely filled and the paste moulded to follow the natural contours of the bone. The Calcibon Inject in the defect should be in direct contact with all surfaces of the defect. Where a load may be applied to the implant it should be used in combination with rigid fixation devices.

#### Step 3:

Secure the surgical site after implanting to prevent motion and any implant displacement. When excess fluid is present at the site, cauterisation, and suction

may be used to reduce bleeding.

### Storage

Store Calcibon Inject between 5°C and 25°C. Direct contact with sunlight or heating systems should be avoided.

### Shelf Life and Disposal

The expiration date is printed on the labelling. **Do not use Calcibon Inject after the expiration date.** Calcibon Inject is environmentally friendly. No special disposal is required. The used applicator should be disposed of as clinical waste.

**Note:** Calcibon Inject is for professional use only. Responsibility for proper selection of patients, for adequate training, for experience in the choice of Calcibon Inject, all aspects of the surgery, and for the choice of post-operative procedures rests entirely with the surgeon.

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Ceramisis Ltd.
Biomet Deutschland GmbH