REF

3020820401-3, 3021180001-3, 3020830401-3,
3021170001-3, 3020840401-3

ZIMMER BIOMET

Refobacin® Plus Bone Cement

Radiopaque bone cement containing gentamicin

Instructions for use
Gebrauchsinformation
Notice d'utilisation
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Käyttöohje
Brugsanvisning
Bruksanvisning
Instruções para uso
Οδηγίες χρήσης
Návod k použití
Kasutusjuhend
Lietošanas instrukcija
Naudojimo instrukcijos
Használati utasítás
Instrukcja użytkowania
Navodila za uporabo
Návod na použitie
Kullanım talimatları
Instrucțiuni de Utilizare
Инструкции за употреба
使用说明
한국어 사용 지침
Инструкция по применению
Інструкція з експлуатації

Reference: 9878300020
IFU updated: 2018/05
Index: 03
First put on market: 2016

EN	Instructions for use	3
DE	Gebrauchsinformation	10
FR	Notice d'utilisation	18
ES	Instrucciones de uso	26
IT	Istruzioni per l'uso	34
NL	Gebruiksaanwijzing	42
FI	Käyttöohje	50
DK	Brugsanvisning	58
SV	Bruksanvisning	65
PT	Instruções para uso	72
GR	Οδηγίες χρήσης	80
CZ	Návod k použití	88
EE	Kasutusjuhend	95
LV	Lietošanas instrukcija	102
LT	Naudojimo instrukcijos	109
HU	Használati utasítás	116
PL	Instrukcja użytkowania	124
SI	Navodila za uporabo	132
SK	Návod na použitie	140
TR	Kullanım talimatları	147
RO	Instrucțiuni de Utilizare	154
BG	Инструкции за употреба	162
CN	使用说明	170
KR	한국어 사용 지침	176
RU	Инструкция по применению	183
UA	Інструкція з експлуатації	191

Refobacin® Plus Bone Cement

EN

Radiopaque bone cement containing gentamicin

IMPORTANT INFORMATION FOR THE OPERATING SURGEON

Before using a product placed on the market by Zimmer or Biomet, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g., product literature, surgical technique). Zimmer or Biomet is not liable for complications that may arise from the use of the device in circumstances outside of Zimmer or Biomet's control including, but not limited to, product selection and deviations from the device's indicated uses or surgical technique.

Intended use and device properties

Refobacin® Plus Bone Cement is a fast setting polymer, containing gentamicin, for use in bone surgery. The gentamicin content protects against infections caused by bacterial invasion of the implant and the surrounding tissue by gentamicin-sensitive strains. Mixing of the two component system, consisting of a powder and a liquid, initially produces a paste, which is used to anchor the prosthesis to the bone. The hardened bone cement allows stable fixation of the prosthesis and transfers stresses produced on movement to the bone via the large interface. Insoluble zirconium dioxide is included in the cement powder as an X-ray contrast medium. The chlorophyll additive serves as optical marking of the bone cement at the site of the operation.

Pharmacological properties

Gentamicin is a parenteral antibiotic which belongs to the class of aminoglycosides. It is a mixture of the structurally very similar homologues gentamicin C1, C1a and C2, C2a and C2b. Gentamicin's mechanism of action is based on disruption of protein biosynthesis on the

bacterial ribosome through interaction with rRNA and subsequent inhibition of translation. This produces a bactericidal effect.

Composition

The cement powder contains:

poly (methyl acrylate, methyl methacrylate)	85 %
zirconium dioxide	12 %
benzoyl peroxide	1 %
gentamicin sulfate	2 %

The liquid monomer contains:

methyl methacrylate	98 %
N,N-dimethyl-p-toluidine (DmpT)	2 %

In addition: hydroquinone, chlorophyll VIII

The content of active gentamicin is given in the following table

Size	Active gentamicin
20	0.28 g
40	0.55 g
60	0.83 g

Presentation and pack sizes

Refobacin® Plus Bone Cement1x20

Original pack consists of:
1 sachet of 22.46 g powder
1 ampoule of 10 ml liquid

Refobacin® Plus Bone Cement2x20

Original pack consists of:
2 sachets of 22.46 g powder each
2 ampoules of 10 ml liquid each

Refobacin® Plus Bone Cement1x40

Original pack consists of:
1 sachet of 44.92 g powder
1 ampoule of 20 ml liquid

Refobacin® Plus Bone Cement2x40

Original pack consists of:
2 sachets of 44.92 g powder each
2 ampoules of 20 ml liquid each

Refobacin® Plus Bone Cement1x60

Original pack consists of:
1 sachet of 67.38 g powder
1 ampoule of 30 ml liquid

Indications

Refobacin® Plus Bone Cementis used for stable anchoring of suitable joint prostheses in primary operations, with added protection against infection, and also in revision operations resulting from aseptic loosening of the prosthesis and infection of the prosthesis by gentamicin sensitive strains.

Contraindications

Refobacin® Plus Bone Cementmust not be used in cases of known hypersensitivity to gentamicin or to other constituents of the bone cement.

Side effects

After preparation of the prosthesis bed or directly after the implantation of the cement and prosthesis, pressure rise in the bone bed may cause a temporary fall in blood pressure.

In addition to hypotension, pulmonary embolism and cardiac arrest with their potentially fatal consequences have been encountered in rare cases. These cardiovascular and respiratory side effects known as the implantation syndrome are mainly caused by infiltration of bone marrow constituents into the venous vascular system.

The site of the prosthesis should therefore be rinsed thoroughly with an isotonic solution (e.g. physiological saline) before implantation.

To minimize the pressure increase in the bone bed during implantation of the cement and the prosthesis, sufficient drainage is recommended.

In the presence of pulmocardiovascular disturbances, the blood loss should be monitored and anaesthesiological measures may be required e.g. in the event of acute respiratory failure.

The constituents of bone cement may lead to local irritations or hypersensitivity reactions in isolated cases.

Application of gentamicin may, in principle, trigger the typical side effects of this antibiotic, which are in particular:

- damage of the acoustic and the vestibular nerve,
- nephrotoxicity,
- neuromuscular blockade (also see interactions),
- in rare cases paresthesia, tetania and myasthenia (due to hypocalcemia, hypomagnesemia and hypokalemia),
- in rare cases allergic reactions (exanthema, urticaria and anaphylactic reactions).

Refobacin® Plus Bone Cementis considered most unlikely to cause gentamicin overdosage, since only low and short-lived serum concentrations are obtained during the first few postoperative hours from the high local gentamicin concentrations.

Adverse Events

When using acrylic bone cements, the following events have been reported:

- displacement of the prosthesis,
- aseptic loosening,
- superficial or deep wound infection,
- hemorrhage and hematoma,
- pain,
- adverse soft tissue reaction,
- periprosthetic fracture,
- implant fracture,
- malalignment,
- wear of the acetabular component,
- lysis,
- dissociation of the liner, pseudarthrosis.

This list may not be complete as it may change over time. For possible further events please refer also to the instruction for use packaged with the prosthesis used in combination with the cement.

Warnings***Intraoperative caution***

Packaging must not show signs that could indicate a defect in the sterility and/or integrity of the medical device. Do not use a device with a damaged packaging.

Preparation of bone cement

Methyl methacrylate is a volatile flammable liquid. The vapor given off during mixing may cause irritation of the eyes and the respiratory passages, as well as more general feeling of ill health such as nausea and headache. Such symptoms can

be minimized through adequate ventilation or by using a closed mixing system.

Monomer (methyl methacrylate) is a lipid solvent. Direct skin contact with the liquid monomer must be avoided as much as possible, as there have been reports of allergic skin reactions (contact dermatitis). It is recommended to wear additional gloves made of polyethylene (PE) during handling of the cement to give adequate protection of the skin.

Application of bone cement

Care must be taken to ensure complete support of all parts of the device embedded in bone cement to prevent stress concentrations which may lead to failure of the implant.

Complete removal of all bone cement debris is critical to minimize or prevent wear of the articulating surfaces.

Implant fracture due to cement failure has been reported.

Polymerization of bone cement is an exothermic process that creates significant heat.

Related to patients

Patients should be warned of the impact of excessive loading that can result if the patient is involved in an occupation that includes substantial walking, running, lifting or excessive muscle loading due to weight that places extreme demands on the joint and can result in device failure or dislocation.

Patient smoking may result in delayed healing, non-healing and/or compromised stability in or around the placement sites.

Pregnancy

No data are available on the administration of Refobacin® Plus Bone Cement during pregnancy. Gentamicin is known to cross the placenta, but limited experience with gentamicin has not suggested an increased risk of malformations. Ototoxicity and nephrotoxicity in the fetus is a potential hazard, but this has not been confirmed clinically. No untoward effects have been described with aminoglycosides in local treatment during pregnancy.

For these reasons, the use of Refobacin® Plus Bone Cement is advised against during pregnancy, unless the benefits for the mother outweigh the potential risk to the child.

Lactation

No data are available on the administration of

Refobacin® Plus Bone Cement during breast-feeding. Gentamicin is excreted in small amounts in human breast milk. Because of enhanced intestinal permeability in neonates, accumulation and ototoxicity cannot be excluded.

For these reasons, the benefits for the mother should outweigh the potential risk to the child before using Refobacin® Plus Bone Cement during lactation.

Use in children

No data are available on use of bone cements in children.

Interactions

Concurrent administration of muscle relaxants and ether may potentiate the neuromuscular blocking properties of gentamicin. However, due to the low serum concentrations this is unlikely to happen.

Incompatibilities

Aqueous (e.g. antibiotic containing) solutions must not be mixed with the bone cement, as this reduces the strength considerably.

Any impurities (e.g. blood, bone debris, glass debris from the ampoule), additional components or any kind of foreign matter should not be mixed with the cement paste.

MR safety

None of the components of Refobacin® Plus Bone Cement include conductive, metallic or magnetic materials. Therefore, in accordance with the definition stated in ASTM-F2503, Refobacin® Plus Bone Cement is classified as "MR Safe" – an item that poses no known hazards in any MR environment.

Dosage

The amount of bone cement required depends on the patient's anatomy and on the implant used. One or more complete units (the contents of one sachet and one ampoule) must always be mixed together in order to ensure the intended handling and mechanical properties of the bone cement. This device has not been tested for more than four 40 g dose units. Even with this dose only low serum levels of antibiotics were determined in clinical studies.

Notes on use

Before using Refobacin® Plus Bone Cement the surgeon must be thoroughly familiar with its

properties, handling and use, and must have read the relevant literature. For special mixing and application techniques (e.g. vacuum mixing, vacuum application, use of a femoral cement restrictor), the relevant instructions of these mixing and application system must be consulted. Before using Refobacin® Plus Bone Cement for the first time, a test-mixing should be carried out to become familiar with the characteristics of Refobacin Plus Bone Cement. The protective aluminium packaging and outer non-sterile polyethylene sachet ("peel-off" package) should be opened, maintaining sterility, by the circulating nurse. The ampoule and powder sachet are then taken out under sterile conditions and placed on a sterile working surface (Figs. 1-4). The ampoule is opened by breaking the neck and the inner sachet is cut open with sterile scissors.

Preparation Method

Refobacin® Plus Bone Cement can be mixed by two different methods:

- Mixing without vacuum,
- Mixing under vacuum.

Mixing the cement under vacuum reduces the porosity of the cement and improves the mechanical stability of the cement.

Mixing without vacuum

The mixing system should be of inert material intended for preparation of PMMA bone cement. Pour the liquid into a sterile mixing bowl and add all the powder.

Stir carefully with a sterile stirring rod or spatula until a homogeneous mass is obtained. The mixture should not be stirred longer than 30 sec, irrespective of the room temperature. The temperature/time diagrams must be consulted (Tab. 1 and Fig. 6).

Vacuum mixing

To obtain a bone cement with lower porosity the use of a vacuum mixing system is recommended (Fig. 5). This method requires an airtight closed system and rapid generation of sufficient vacuum (approx. 200 mbar absolute pressure) in the mixing cartridge. The stirring time for vacuum mixing is the same as for mixing without vacuum (30 sec), (Tab. 2 and Fig. 7). Method of mixing is given in the instructions for the system used and must be consulted.

Application procedure

After careful preparation of the bone bed, Refobacin® Plus Bone Cement can be applied manually or using a cement syringe or some


Fig. 1


Fig. 2


Fig. 3

Fig. 4


Fig. 5

alternative application system. Method of application is given in the instructions for the system used and surgical technique for the prosthesis being implanted and both must be consulted. The working time and the rate of polymerization are strongly dependent on ambient temperature and on temperature of the components. The hardening time is shorter at high temperatures and longer at low temperatures. Viscosity increases as polymerisation progresses, i.e. as the cement hardens.

Manual Application timing

Depending on the ambient temperature and the temperature of the components, the following schedule applies (Tab. 1). Fig. 6 shows the duration of sticky and dough phases as a function of temperature. The cement and prosthesis must be applied before the end of dough phase. Applying the cement after the end of the dough phase may cause uneven and/ or inadequate cement-bone bonding, thereby resulting in possible early loosening of the implant (see graphs). After positioning of the implant any movement must be avoided, to ensure the anchoring of the prosthesis.

Application timing of vacuum-mixed cement

Vacuum-mixed cement is usually applied with a syringe. Depending on the ambient temperature and the temperature of the components the following schedules apply: Tab. 2 and Fig. 7. The instructions of the mixing system used must be consulted.

Use in joint surgery

An appropriate up to date cementing technique must be used with Refobacin® Plus Bone Cement to limit side effects (e.g. embolism) and to achieve a stable and long-lasting anchorage of the prosthesis. A pre-requisite for this is careful preparation of the prosthesis site by thorough rinsing (e.g. with physiological saline) before application of the cement. Sufficient drainage is recommended to avoid any pressure build-up in the bone bed during implantation. The cementing technique should be adapted to the shape of the implanted site.

Storage

Do not store above 25°C.
Keep away from sunlight.
Keep dry.

Shelf life/sterility

The expiry date is printed on the outer carton, on the protective aluminum packaging and on the inner sachet. Refobacin® Plus Bone Cement must not be used after the expiry date. The contents of unused but opened or damaged packs must not be re-sterilized and are to be discarded.

This product is for single use only and has not been designed or tested for re-use. Do not attempt to clean or re-sterilize the product due to the risk of cross infection and/ or possible alterations in the product performance.

The liquid monomer has been sterilized by aseptic processing. The powder and the outside of the monomer ampoule have been sterilized with ethylene oxide. An expiry date is printed on the packaging of the liquid monomer.

Refobacin® Plus Bone Cement must not be re-sterilized.

If the polymer powder shows a yellow discoloration, it must not be used.

Disposal

The contents of unused but opened packaging must not be re-sterilized and are to be discarded. Disposal of monomer contents/container is to be performed in accordance with local regulation. Any disposal of powder, and liquid, cured bone cement or mixing system should be treated as biological waste.

PATIENT WARNINGS

Postoperative care is important. The patient must be warned of the limitations of the reconstruction and the need for protection of the implants from full weight bearing or load bearing until adequate fixation and healing has occurred. The patient is to be cautioned to govern activities, protecting the joint replacement from unreasonable stress conditions. Excessive activity, failure to control body weight and trauma affecting the joint replacement have been associated with premature failure of the reconstruction by loosening, fracture and/ or wear of the implants. Loosening of the components can result in increased production of wear particles, as well as accelerating damage to bone, making successful revision surgery more difficult.

The patient is to be advised of the importance of postoperative follow-up examination and that excessive activity or trauma can lead to failure.

The patient is to be made aware of and warned in advance of the surgical risks and of possible adverse effects. The patient must be warned that the device does not replace normal healthy bone, and that the implant can break or be damaged as a result of excessive load bearing or trauma. The patient must be warned to inform any other medical practitioner who may treat him in the future of the presence of the implant.

Responsible manufacturer

Biomet France
Plateau de Lautagne
26000 Valence
France
Phone: 0033 4 75 75 91 00

CE 0123

Information status 2018/05

REFOBACIN® is a registered trademark owned by and licensed from Merck KGaA.

Tab. 1: Manual application of REFOBACIN® PLUS BONE CEMENT

Ambient and component temperature	18°C	20°C	23°C
End of mixing at ambient temperature	0'30''	0'30''	0'30''
End of sticky phase	3'00''	2'30''	1'45''
End of dough phase	7'00''	5'45''	4'45''
Hardening	13'00''	11'15''	9'15''

Tab. 2: Application of REFOBACIN® PLUS BONE CEMENT mixed under vacuum

Ambient and component temperature	18°C	20°C	23°C
End of mixing at ambient temperature	0'30"	0'30"	0'30"
End of sticky phase	2'00"	1'45"	1'30"
End of dough phase	5'30"	4'45"	3'45"
Hardening	12'00"	10'15"	8'30"

Fig. 6: Manual application of REFOBACIN® PLUS BONE CEMENT



Fig. 7: Application of REFOBACIN® PLUS BONE CEMENT mixed under vacuum

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